DV51 STANDARD CPAP SERIES

EN DeVilbiss® IntelliPAP® Positive Airway Pressure Device

CAUTION—USA Federal law restricts this device to sale by, or on the order of a physician.

Assembled in USA of US & Imported Parts. Instruction Guide.

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








SYMBOL DEFINITIONS

	Class II electrical protection-double insulated		Type BF equipment-applied part
	Data port input/output		Standby-turns blower on/off
	Next item on LCD display		Increase value on LCD display
	Previous item on LCD display		Decrease value on LCD display
	Attention, consult accompanying documents		Heat
	DC input: 12 VDC, center pin positive		100-240 VAC input 50/60 Hz
	Locked		Unlocked
	IPX1 Drip-proof, vertical		Delay-starts the delay cycle
	This device contains electrical and/or electronic equipment that must be recycled per EC Directive 2002/96/ EC-Waste Electrical & Electronic Equipment		

IMPORTANT SAFEGUARDS

READ ALL INSTRUCTIONS BEFORE USING THIS DEVICE.

SAVE THESE INSTRUCTIONS.

DANGER

- **Electric Shock Hazard** – Do not use while bathing.
- **Electric Shock Hazard** – Do not immerse this device into water or any other liquid.
- **Electric Shock Hazard** – Do not attempt to open or remove the cabinet. There are no user-serviceable internal components. If service is required, return the product to your home care provider. Opening or tampering with the product will void the warranty.

WARNING

- The DeVilbiss IntelliPAP® should be used only with masks recommended by DeVilbiss, your physician or respiratory therapist.
- To avoid rebreathing of exhaled air, do not use a CPAP mask unless the device is turned on and providing a supply of air. Venting in the mask should never be blocked. When the device is turned on and providing a fresh supply of air, exhaled air is flushed out of the mask vent. However, when the device is not operating, exhaled air may be rebreathed. Rebreathing of exhaled air for longer than several minutes can in some circumstances lead to suffocation. This warning applies to most CPAP devices.
- The DeVilbiss IntelliPAP is not a life support device and may stop operating with certain device faults or with a power failure. It is intended to be used on spontaneously breathing individuals weighing 66 lbs/30 Kg or greater.
- To avoid electric shock, always unplug power cord from wall outlet power source when performing cleaning.
- Use only accessories recommended by DeVilbiss.
- Always remove the water chamber from the optional humidifier cradle before filling.
- If water has spilled onto the humidifier cradle or it has been submerged into water, unplug power cord from power source immediately. Allow device to dry completely before use.
- Never block air openings of the humidifier cradle or chamber. Do not insert objects into any openings or tubes.
- For proper operation, place system on a flat, sturdy, level surface.
- The optional humidifier is intended for single-patient use only.
- If the system is placed on the floor, do not position it where it may be stepped on when getting out of bed.

CAUTION

- The circular data port connector located on the back of the IntelliPAP is used to attach accessories to the device. The connector must only be used with accessories approved for use by DeVilbiss. Do not attempt to attach any other device to this connector as it may damage the CPAP or the accessory device.

- Never rinse or place the device in water. Never allow liquids to get into or around any of the ports, switches or air filter; doing so will result in device damage. If this occurs, discontinue use and remove the power cord from the power source. Allow the device to completely dry before use.
- Do not place the IntelliPAP where it can be knocked onto the floor or where the power cord may create a trip hazard.
- Only the DeVilbiss DV5 series Heated Humidifier system is recommended for use with the IntelliPAP. Other humidifier systems may prevent the device from detecting snoring and may cause inappropriate pressure levels in the mask.
- Empty and dry humidifier water chamber before transporting.

INTRODUCTION

Intended Use

The DeVilbiss® IntelliPAP® Model DV5 Series is intended for use in treating obstructive sleep apnea in patients 66 lbs/30 kg and above.

Indications For Use (Optional Heated Humidifier)

Use on the advice and prescription of a licensed physician to help relieve the symptoms of dryness of the throat, nasal passages and the mouth, which are common with positive airway pressure therapy. This is especially true in dry climates and during the cold season when humidity in the air is typically lower than at other times.

Contraindication

Do not use IntelliPAP system if your upper airway has been bypassed.

NOTE—*Symptoms of dryness of the throat, nasal passages, and the mouth are common with positive airway pressure therapy. The IntelliPAP features an optional humidifier system to help minimize these effects.*

Key Features

DeVilbiss IntelliPAP Device (Figures A and B)

1. Keypad/LCD Display (see Figure C)
2. Air Supply Port on Back
3. Air Supply Port on Bottom (for optional humidifier)
4. Air Supply Port Plug
5. Heater Power Connector (for optional humidifier)
6. AC Power Connector
7. DC Power Connector
8. Data Port
9. Power Cord (appropriate to your wall outlet) (not shown)
10. Air Inlet Filter Opening
11. Connector Cover (not shown)

Keypad (Figure C)

1. On/Off
2. Previous Item
3. Next Item
4. Delay
5. Decrease Value
6. Increase Value
7. Heater Power LED (for optional humidifier)

Ask your equipment provider for information about other DeVilbiss masks, equipment, and accessories.

Humidifier Cradle (Figure D)

1. Heater Plate
2. Heater Power Connector
3. Storage Compartment for Connector Cover (on bottom)
4. Air Supply Port Inlet
5. Air Supply Port Outlet
6. Flow Generator Release Button
7. Flow Generator Locking Tabs

Water Chamber (Figure E)

1. Chamber Lid
2. Chamber Base
3. Heat Transfer Plate (on bottom)
4. Sealing Gasket
5. Water Level Indicators (front and sides)
6. Chamber Release Latch
7. Chamber Disassembly Lever

SYSTEM ASSEMBLY

Without Heated Humidification

1. Locate the air supply port and heater power connector on the bottom of the IntelliPAP device. Ensure that each opening has the appropriate soft cover securely attached.
2. Place the IntelliPAP on a stable surface such as a nightstand or table. If you wish, you may also place the device on the floor at the bedside, being careful to place the unit where it won't be kicked or stepped on. Ensure that the air inlet in the back of the device is not blocked by anything such as curtains or bedding.


NOTE—*Never place the IntelliPAP system on a soft surface such as a bed or couch during operation.*

With Heated Humidification

CAUTION–The height of the IntelliPAP must be lower than the mask when using a humidifier to prevent water from getting into mask.

1. If your system includes a humidifier, please follow steps in Figures 1-7.
 - a. Remove the air supply tubing from the back of the device.
 - b. Locate the air supply port at the bottom of the device. Remove the port plug and insert it into the air supply port in the back of the device.
 - c. Locate the heater power connector on the bottom of the device. Remove the connector cover and insert it into the storage opening in the bottom of the humidifier cradle.

KEYPAD AND DISPLAY

Basic CPAP Controls		
Name	Symbol	Function
ON/OFF		The ON/OFF button starts and stops the flow of air from the device.
HEATER		If your IntelliPAP has a heater/ humidifier attached, these buttons adjust the heater setting. Increase the setting for more humidity and decrease the setting for less humidity.
DELAY		The Delay button activates the comfort delay feature and only functions if this feature has been enabled. To understand how to set the parameters of this feature, please refer to the Using the Comfort Delay Feature found on page EN - 6.

NOTE– The left and right arrow keys on the control panel provide access to menus for specific settings. Ask your homecare provider for more information.

NOTE– Advanced directions (ADV-DV5X) for utilizing device features are available for download on the Product Literature section of our website or can be ordered by calling Customer Service.

OPERATION

1. Fully insert the power cord into the back of the device. Plug the other end of the power cord into the AC wall outlet.
2. Attach air supply tubing to air supply port in the back of the device.
3. Prepare the mask (supplied by your equipment provider) by referring to the mask's instruction guide.
4. Connect the air supply tubing to the mask.

NOTE–The air supply tubing supplied with the device is specially designed to prevent kinking and allow you to receive optimal airflow. Use only the 22 mm diameter, smooth-bore, kink-resistant tubing supplied with this device or an air supply tubing recommended by your equipment provider.

Start-up

Put on your mask and press the ON/OFF button on the keypad OR put on your CPAP mask and take a few breaths to turn on airflow automatically. The pressure in your mask will rise to the set pressure in about 10 seconds.

NOTE–If the device does not automatically start after taking 2-3 breaths, this feature may have been disabled. If Auto ON/OFF is disabled, you must manually press the ON/OFF button to start and stop the device.

The keypad display shows the actual mask pressure while the device is in use. Because this is the actual pressure, the device pressure reading may vary slightly as you breathe.

If your DeVilbiss IntelliPAP system includes a humidifier, its heater controls and status are displayed.

NOTE– If the display on your device is not similar to those listed above, refer to the section titled "Patient Messages" on page EN-8.

Using the Comfort Delay Feature

Your IntelliPAP has a comfort delay feature to help make it easier for you to fall asleep. The airflow will start at a low pressure and then gradually ramp-up to the prescribed setting over the last few minutes of the specified delay time.

NOTE–When the delay comfort feature is activated, a message will flash on the display periodically, indicating the amount of time remaining in the delay.

You can adjust the amount of time spent in comfort delay by using the advanced device controls. The Delay Time can be adjusted from 0 to 45 minutes in 5-minute increments.

NOTE–If you are unable to change the delay time, your provider has locked this setting.

Press the delay button to activate the comfort delay feature. The airflow pressure will drop to the prescribed comfort delay pressure.

In the final 10 minutes of the delay period, the pressure will gently rise to the pressure prescribed for your therapy. You may repeat the comfort delay by pressing the delay button again if you have not fallen asleep during the first use.

Controlling the Humidity (optional)

The heater setting can be adjusted within a range of 1 to 10. The optimum temperature setting is dependent on the ambient temperature and humidity of your room. If your room temperature is low and/or the relative humidity is high, then a lower temperature setting may be necessary to prevent excess condensation in the air supply tubing. If the room temperature is high and/or the relative humidity is low, then a higher temperature setting may be necessary for optimum comfort.

The heater setting of 1 results in a heater plate temperature of approximately 84°F (29°C). The heater setting of 10 results in a heater plate temperature of approximately 149°F (65°C).

NOTE—Always allow at least 10 minutes of cooling time after turning the device off before removing the water chamber from the humidifier cradle.

WARNING

Do not attempt to fill the humidifier while it is attached to the humidifier cradle. Damage to the humidifier cradle may occur.

CAUTION—Use only distilled water that is at room temperature. Do not add any medication or other additives to the water.

CAUTION—For proper operation, ensure the flow generator is OFF before inserting humidifier chamber.

WARNING

Never touch the heater plate on the humidifier cradle. Never touch the heat transfer plate on the bottom of the water chamber. These plates can reach temperatures as high as 149°F (65°C) during operation.

Do not operate the heater if the water chamber is empty. The heater plate may be turned off using the keypad heater control when being used without water.




Shut Down

Press the ON/OFF button to turn the device off. You may also simply remove your mask. After a few seconds the device will display the message “Mask Off. Please check mask fit.” If you do not put the mask on again the device will automatically turn off in about 20 seconds. Pressing any button will remove the message from the display.

NOTE—If the DeVilbiss IntelliPAP does not turn off automatically, the Auto-OFF feature may have been disabled. In this case you must manually turn the device off.

Enable Menu List

The Enable menu is available in all modes and is a sub-menu that allows you to control various functions of your IntelliPAP. It is accessed by pressing the down button while “Enable Menu” is shown on the Display.

Display Example	Enable Menu Item Description
	This option controls the Auto-OFF feature, which automatically stops the flow of air after taking off the mask. The device will require at least 30 seconds to turn OFF depending on the kind of mask used. When Auto-OFF is enabled, Auto-ON is also enabled. If the unlocked symbol is shown, you are able to adjust this setting. If the locked symbol is shown, this setting can only be adjusted by your provider.
	This option controls the Auto-ON feature, which automatically starts the flow of air after breathing once or twice into the mask. This is always enabled when Auto-OFF is enabled. If the unlocked symbol is shown, you are able to adjust this setting. If the locked symbol is shown, this setting can only be adjusted by your provider.
	Low backlight Enable will keep the backlight dim during operation. Low backlight Disabled will turn the backlight OFF during operation.
	Pressing the up key exits the Enable menu.

SmartCode and Adherence Score Information

Your homecare provider may contact you to retrieve SmartCode and Adherence Score information. To display the SmartCode, press the LEFT arrow key once. Use the UP arrow or DOWN arrow keys to change the SmartCode reporting period. Press the LEFT arrow key again to display the Adherence Score.

PATIENT MESSAGES

Your IntelliPAP will alert you of issues that may require some action on your part in order to make your therapy more effective. You can press any key or button to clear the message from the display, but the IntelliPAP will operate normally even if the message is not cleared. Two kinds of messages are displayed: Reminders and Notifications.

Reminders

Reminders alert you when parts of your system need to be replaced. Some components of your system wear over time and, if not replaced, may compromise your therapy. If you see a message on the display that is not listed below, you may have the optional DeVilbiss SmartLink™ Module attached to your device, which provides additional messages. Refer to the documentation that was provided with the SmartLink Module.

Clean Filter – The filter should be checked every 10 days for signs of dirt or wear so that it can be cleaned as needed. This message helps to remind you to check it regularly.

Notifications

Notifications identify conditions in your IntelliPAP that may require action by you or your equipment provider in order to maintain a high level of therapy. If you see a message on the display that is not listed below, you may have the optional DeVilbiss SmartLink Module attached to your device, which provides additional messages. Refer to the documentation that was provided with the SmartLink Module.

Delay Running xx Minutes Left – While a comfort delay session is active, this message is flashed on the display every 5 seconds to let you know how much time is left in the delay.

Mask Off – This message appears due to a poor mask fit or a removed mask. Check for leak around the mask seal and make adjustments as necessary according to the mask manufacturer's instructions. This message will be displayed until the high airflow problem is corrected. If the high airflow condition persists for about 20 seconds, the device will automatically turn off if Auto-OFF is enabled.

Device Fault – If a device fault message is displayed, refer to Troubleshooting for instructions.

INTELLIPAP TRAVEL INFORMATION

International Power Changes

Your DeVilbiss IntelliPAP is equipped with a universal power supply. This means that it is automatically capable of accepting line voltages of 100 to 240V~ 50/60Hz. You do not need to make any adjustments to the device if you are traveling to another country where the line voltage is different than your domestic line voltage. Simply contact your equipment provider for the correct style power cord for the region in which you will be traveling (refer to Accessories/Replacement Items for the correct part number).

Altitude

Your IntelliPAP automatically compensates for changes in altitude between sea level and 9000 ft (2750 m). You do not need to make adjustments to the device for changes in altitude.

INTELLIPAP DC OPERATION

Battery Power

Your IntelliPAP can be operated from a 12V DC power source in a motor home, boat or other recreational vehicle. The power source can be used with a 12V DC power cord or a 12 V DC stand-alone battery. Both parts are available from your equipment provider (refer to Accessories/Replacement Items for the correct part numbers). Both accessory cables are required for stand-alone battery operation. The stand-alone battery connection uses clamps that attach to the terminal connections of the battery.

NOTE–The heater for the IntelliPAP's humidification system will not operate directly from a 12V DC power source. If you need humidification, non-heated humidification is still available. Alternately, instead of connecting 12V DC power directly to the IntelliPAP, you can pass the DC power through an inverter that will deliver AC power to the system. This will allow the heater to function. The power rating of the inverter must be at least 200 watts @ 110V~ / 400 watts @ 220V~.

NOTE—The IntelliPAP will automatically switch power sources between AC and 12V DC if both sources are available. If your unit is connected to both AC and 12V DC power sources, it will operate on AC power as long as AC power is present, and automatically switch to 12V DC power if AC power is lost. When AC power is restored, it will return to the AC source. You do not need to reset or adjust any settings when the power source switches, but remember that your heater will not operate on 12V DC.

SUPPLEMENTAL OXYGEN

WARNING-OXYGEN USE

- Oxygen supports combustion. To avoid possible physical injury, do not smoke while using this device with supplemental oxygen. Do not use this device near hot objects, volatile substances or sources of open flames.
- Always turn on the IntelliPAP before turning on the oxygen source. Turn off the oxygen source before turning off the IntelliPAP.
- This device is not suitable for use in the presence of a flammable anesthetic mixture.
- At a fixed flow rate of supplemental oxygen, the inhaled oxygen concentration will vary depending on the pressure setting, patient breathing pattern, mask selection, and leak rate. This warning applies to most types of CPAP devices.
- Never allow the oxygen source to run continuously while connected to the IntelliPAP if the device is not in use. If the IntelliPAP is not being used, turn off the oxygen flow.

CAUTION

- Oxygen is a prescribed gas and should only be administered under the supervision of a physician.
- The flow setting of the oxygen source must be specified by a physician.

Supplemental oxygen is sometimes prescribed in addition to CPAP for treatment of Obstructive Sleep Apnea. If your physician prescribed oxygen for your sleep therapy, you can add the oxygen in one of two ways:

1. Obtain the optional oxygen adapter (part #7353D-601) from your homecare provider and connect it to the air supply port on the IntelliPAP. Connect one end of the air supply tubing to the open end of the adapter and the other end of the air supply tubing to the mask.
2. If your mask is equipped with an oxygen port, your physician may request that you apply oxygen directly to your mask.

ACCESSORIES/REPLACEMENT ITEMS

Replacement Items			
Air-inlet filter (4/pk)	DV51D-602	Heater connection cover	DV51D-605
Fine particle filter (4/pk)	DV51D-603	Oxygen adapter	7353D-601
Air supply tubing (6')	DV51D-629	Carrying case	DV51D-610
Air supply plug	DV51D-604	Heated Humidifier	DV5HH
Humidifier Chamber	DV5C	Chamber Sealing Gasket	DV5C-614
Power Items			
AC power cord (USA)	DV51D-606	AC power cord (Australia)	DV51D-609
AC power cord (UK)	DV51D-608	DC power cord	DV51D-619
AC power cord (EU)	DV51D-607	DC battery clamp-on adapter (requires DC Power Cord DV51D-619)	DV51D-696
DeVilbiss Masks (USA) Recommended for use with IntelliPAP			
FlexSet	9354D	Serenity	9352D
FlexSet Shallow	9354S	Serenity Shallow	9352S
FlexSet Gel	9354G	Serenity Gel	9352G
FlexSet Gel Shallow	9354GS	Serenity Gel Shallow	9352GS
EasyFit Nasal	Available in Silicone (S, M, L, XL) or Gel (S, M, L). Contact your equipment provider		
EasyFit Full Face	Available in Silicone (S, M, L, XL) or Gel (S, M, L). Contact your equipment provider		

MAINTENANCE

DANGER

Electric shock hazard-Do not attempt to open or remove the cabinet; there are no user-serviceable components inside. If service is required, contact your equipment provider for instructions on obtaining service. Opening or attempting to service your device will void the warranty.

WARNING

Disconnect supply before servicing.

Pressure Accuracy—The IntelliPAP is designed to provide accurate pressure performance from the factory without additional calibration. No routine service is required provided the device is used in accordance with the manufacturer's directions.

NOTE—Some countries require periodic servicing and calibration of this kind of medical device. Contact your equipment provider for further information.

Standard Air-Inlet Filter—Check the filter every 10 days and clean it if necessary. Replace the filter every 6 months or sooner if it is damaged. The standard filter is designed to filter particles larger than 3.0 microns in size. The IntelliPAP MUST have the standard filter in place during operation.

Optional Fine Particle Filter—Check the filter every 10 days and replace it if it is dirty or damaged; otherwise replace the filter every 30 days. The optional fine particle filter is designed to filter particles as small as 0.3 microns in size.

CLEANING

WARNING

To avoid electric shock, always unplug the power cord from the wall outlet power source before cleaning.

CAUTION—Never rinse or place the device in water. Never allow liquids to get into or around any of the ports, switches or air filter; doing so will result in product damage. If this occurs, do not use the device. Contact your equipment provider for service.

Device Cabinet

1. Unplug the IntelliPAP from the power source. Wipe the cabinet with a clean, damp cloth every few days to keep the cabinet dust free.
2. Allow the device to dry completely before plugging it into a power source.

Filters

Air-Inlet Filter

1. The air-inlet filter should be checked every 10 days and cleaned as necessary. The filter is located in the rear of the device.

CAUTION—Proper filter function is important for the operation of the device and to protect the device from damage.

2. To clean, remove the dark outer foam filter from the frame.
3. Wash the filter in a solution of warm water and dishwashing detergent, and rinse with water. Allow the filter to dry completely. If the filter (part #DV51D-602) is damaged, contact your equipment provider for a replacement.
4. Ensure the filter is completely dry before installation and use.

Fine Particle Filter (Optional)

1. If you have an optional fine particle filter (part #DV51D-603), inspect this filter and replace it if it shows signs of dirt or damage. Do not wash.
2. Install the fine particle filter behind the standard filter.

NOTE—If the fine particle filter is not installed first according to the directions, the life of the filter will be reduced requiring more frequent replacement.

Optional Heated Humidifier

Water Chamber - Humidifier water chamber should be cleaned daily.

1. Press on the chamber release latch and pull out the water chamber from the cradle (Fig 4).
2. Separate the chamber halves by pressing forward on the chamber disassembly lever and lifting the lid of the chamber (Fig 9).
3. Remove the sealing gasket from the chamber lid (Fig 10).
4. Wash the parts in a solution of warm water and mild dishwashing detergent.

CAUTION—Hand wash only. Do not use any solutions that contain bleach, chlorine, alcohol, scented oils or other aromatic-based ingredients to clean the humidifier components or tubing. These solutions may reduce the life of the product.

5. Allow the parts to dry before reassembly.

CAUTION—Inspect the gasket seal before each use to ensure that it is properly seated in the groove around the entire chamber.

Humidifier Cradle - The humidifier cradle should be cleaned as needed.

1. Turn off the heater and unplug from power source. Allow system to cool for at least 10 minutes.
2. Detach flow generator from the cradle; first remove water chamber then fully depress the flow generator release button to separate (Fig 4).

NOTE—Flow generator release button must be pressed completely to separate.

3. Use a soft, damp cloth to clean the surfaces of unit.

Tubing

The air supply tubing should be cleaned daily. Remove the tubing from the device and mask. Use a mild detergent and water to clean the inside of the tubing. Rinse the tubing and allow to air dry.

Mask and Headgear

Clean the mask and headgear according to the manufacturer's instructions.

TROUBLESHOOTING

WARNING

Electric Shock Hazard - Do not remove the device cabinet; there are no user-serviceable internal components. The cabinet should only be removed by a qualified DeVilbiss provider.

Issue	Possible Cause	Remedy
Nothing shows on the display.	1. Device is not plugged in or the power cord is not fully inserted.	1a. Verify that power cord is firmly connected to the IntelliPAP and the power outlet. 1b. If you are using a DC power source, make sure the cable connections are secure. Verify that battery is charged.
	2. Optional Humidifier-The flow generator is not fully seated on humidifier cradle.	2. Refer to assembly instructions (Fig. 1-7).
	3. There is no outlet power.	3. Locate a suitable functioning power source.
The IntelliPAP does not start when breathing into the mask.	1. The Auto-ON feature has been disabled.	1. Use the ON/OFF button to start and stop device.
	2. There is no power to the device.	2. Verify that the power cord is firmly connected to the device and the power outlet.
	3. Breath is not deep enough for the Auto-ON feature to detect.	3. Breathe deeply in and out to start the device.
	4. You are using a full face mask that has an anti-asphyxia valve.	4. Auto-ON may not work because your exhaled breath escapes through the open valve. Use the ON/OFF button to start and stop the device.
	5. The air supply port plug is missing or not fully inserted.	5. Make sure the air supply port plug is fully inserted into the unit.
	6. The humidifier chamber is not fully engaged into the cradle or is missing.	6. Slide the humidifier chamber fully into the slot. Verify the latching lever snaps into place.
The airflow has stopped unexpectedly during use or reports mask leak.	1. The Auto-OFF feature has detected large airflow because of loosely fitting mask.	1. Ensure you have a good mask fit; adjust mask and headgear.
	2. During use, your mouth opens and you begin to mouth-breathe.	2. Contact equipment provider for a chin strap or different mask to prevent mouth-breathing.
	3. The optional humidifier chamber is not fully engaged or is missing.	3. Slide the humidifier chamber fully into cradle. Verify the release latch clicks into place.
	4. The air supply port plug is missing from the back of the CPAP or is not fully inserted.	4. Ensure the air supply port plug is fully inserted into the back of the CPAP.

Issue	Possible Cause	Remedy
The display shows a device fault.	An error has occurred in the device and requires servicing.	Contact your equipment provider for service.
Apnea symptoms have recurred.	1. Air filter may be dirty.	1. Clean or replace the air filter and relocate device away from drapes or other dusty surfaces.
	2. Apnea condition has changed.	2. Contact your physician or equipment provider.
Skin becomes irritated where mask contacts face.	1. Headgear is too tight or improperly adjusted.	1. Loosen headgear to reduce contact pressure on face.
	2. Your mask may not be sized properly or is not the most appropriate shape for you.	2. Contact physician or equipment provider.
Dryness of throat or nose.	Inadequate humidity.	1a. Add a DeVilbiss DV5 heated humidifier. 1b. Increase the heater setting on the humidifier.
Water condensation collecting in the hose causing a gurgling sound.	1. Humidification is too high.	1. Reduce humidifier heater setting.
	2. Room temperature fluctuates from higher to lower levels through the night.	2. Increase room temperature.
Air from flow generator seems to be too warm.	1. Air filters are dirty.	1. Clean filters.
	2. Air inlet port is blocked.	2. Unblock air inlet.
	3. Room temperature is too high.	3. Lower room temperature.
	4. Device is located near a heat source.	4. Locate device away from heat source.
	5. Humidifier (optional) heater setting is too high.	5. Lower heater setting.
Nasal, sinus or ear pain, runny nose	You may have a reaction to air flow pressure.	Discontinue use and contact physician.
Heater setting is not visible on display	1. The flow generator is not fully seated onto the humidifier cradle.	1. Ensure proper contact is made (Fig. 3).
	2. Your CPAP device is not equipped with the optional humidifier.	

SPECIFICATIONS

CPAP

Size.....	4.2"H x 6.5"W x 6.9"D (10.7cm x 16.5cm x 17.5cm)
Weight.....	1.22 kg (2.7 lbs.)
Electrical Requirements AC.....	100-240V~, 50/60 Hz
Electrical Requirements DC.....	10.4 - 15 VDC, 5 Amps
Maximum Power Consumption.....	.65 watts max from AC power source (flow generator only)
Pressure Range.....	3-20 cmH ₂ O
Operating Temperature Range.....	41°F to 104°F (5°C to 40°C)
Operating Humidity Range.....	0% to 95% RH non-condensing
Operating Atmospheric Conditions.....	Sea level to 9,000 feet (2743 m)
Storage & Transportation Temperature Range.....	-40°F to +158°F (-40°C to 70°C)
Storage & Transportation Humidity Range.....	0% to 95% RH non-condensing
Maximum Limited Pressure.....	20 cmH ₂ O under normal use
Sound Pressure Level (tested per ISO 17510).....	26 dBA
Sound Power Level.....	38 dBA

Humidifier

Humidity Output (in the operating flow range):.....	≥10 mgH ₂ O/l air
Size:.....	2.6" H x 6.3" W x 8.4" D (6.6 cm x 16.0 cm x 21.3 cm)
Weight:.....	1.75 lbs. (0.794 Kg)

Power Consumption: 85 Watts
 DV5HH Electrical Requirements 100-240VAC, 0.95 A
 Heater Plate Thermal Cutout: 221°F (105°C)
 Power/Temperature Control: Setting of 1 (minimum) to 10 (maximum)
 Heater Plate Temperature: 84°F to 149°F (approx. 29°C to approx. 65°C)
 Operating Atmospheric Pressure Range: 70.0 – 106.0 kPa

Filter Specifications

Standard Filter > 3.0 micron particles
 Optional Fine Particle Filter > 0.3 micron particles

Dynamic Short-Term Pressure Accuracy with and without DV5HH Humidifier (Per ISO 17510-1:2007)

Waveform	Volume [mL]	Breath Rate [min-1]	Set Pressure [cmH ₂ O]	Pressure Accuracy [Pk-Pk cmH ₂ O]
Sinusoidal cycle with I:E = 1:1	500	10, 15, 20	3.0	0.5
			7.0	0.5
			12.0	1.0
			16.0	1.0
			20.0	1.0

Maximum Flow-Rate (per ISO 17510-1:2007)

Test Pressures					
	3.0 cmH ₂ O	7.0 cmH ₂ O	12.0 cmH ₂ O	16.0 cmH ₂ O	20.0 cmH ₂ O
Measured pressure at the patient connection port (cmH ₂ O)	2.1	6.0	11.0	15.0	19.0
Average flow at the patient connection port (l/min)	70.9	81.5	86.3	90.1	93.1

Pressure Swings (cmH₂O peak to peak at 500ml tidal volume, sine wave profile)

Breaths per minute	10	15	20
@ 6.5 cmH ₂ O	0.5	0.5	0.5
@10 cmH ₂ O	1.0	1.0	1.0
@13 cmH ₂ O	1.0	1.0	1.0
@20 cmH ₂ O	1.0	1.0	1.0

Long Term Static Pressure Accuracy +/- 0.5 cmH₂O

Additional Specifications

Warranty Please contact homecare provider
 Equipment classification with respect to protection from electric shock Class II
 Degree of protection from electric shock Type BF Applied Part
 Degree of protection against ingress of liquids IPXI Drip-proof vertical
 Mode of operation Continuous
 Equipment not suitable for use in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide

DEVILBISS GUIDANCE AND MANUFACTURER'S DECLARATION

WARNING

Medical Electrical Equipment needs special precautions regarding EMC and needs to be installed and put into service according to the Electromagnetic Compatibility [EMC] information provided in the accompanying documents.

Portable and Mobile RF Communications Equipment can affect Medical Electrical Equipment.

The equipment or system should not be used adjacent to or stacked with other equipment and that if adjacent or stacked use is necessary, the equipment or system should be observed to verify normal operation in the configuration in which it will be used.

NOTE– The EMC tables and other guidelines provide information to the customer or user that is essential in determining the suitability of the Equipment or System for the Electromagnetic Environment of use, and in managing the Electromagnetic Environment of use to permit the Equipment or System to perform its intended use without disturbing other Equipment and Systems or non-medical electrical equipment.

Guidance and Manufacturer's Declaration – Emissions All Equipment and Systems			
This device is intended for use in the electromagnetic environment specified below. The customer or user of this device should assure that it is used in such an environment.			
Emissions Test	Compliance	Electromagnetic Enforcement – Guidance	
RF Emissions CISPR 11	Group 1	This device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.	
RF Emissions CISPR 11	Class B Radiated and Conducted Emissions	This device is suitable for use in all establishments including domestic, and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.	
Harmonics IEC 61000-3-2	Class A		
Flicker IEC 61000-3-3	Complies		
Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment - Guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	±6kV contact ±8kV air	±6kV contact ±8kV air	Floors should be wood, concrete or ceramic tile. If floors are synthetic, the relative humidity should be at least 30%.
Electrical Fast Transient/burst IEC 61000-4-4	±2kV on AC Mains	±2kV on AC Mains	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1kV Differential ±2kV Common	±1kV Differential ±2kV Common	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	>95% Dip for 0.5 Cycle 60% Dip for 5 Cycles 30% Dip for 25 Cycles >95% Dip for 5 Seconds	>95% Dip for 0.5 Cycle 60% Dip for 5 Cycles 30% Dip for 25 Cycles >95% Dip for 5 Seconds	Mains power quality should be that of a typical commercial or hospital environment. If the user of this device requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible power supply or battery.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment - Guidance
Power Frequency 50/60Hz Magnetic Field IEC 61000-4-8	3A/m	3A/m	Power frequency magnetic fields should be that of a typical location in a typical commercial or hospital environment.
Conducted RF IEC 61000-4-6	3 Vrms from 150 kHz to 80 MHz	V1 = 3 Vrms	Portable and mobile RF communications equipment should be separated from the device by no less than the recommended separation distances calculated/listed below: $D=(3.5V1)\sqrt{P}$
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	E1 = 3V/m	$D=(3.5/E1)\sqrt{P}$ 80 to 800 MHz $D=(7/E1)\sqrt{P}$ 800 MHz to 2.5 GHz Where P is the maximum power rating in watts and D is the recommended separation distance in meters. Field strengths from fixed transmitters, as determined by an electromagnetic site survey, should be less than the compliance levels (V1 and E1). Interference may occur in the vicinity of equipment containing a transmitter.

For transmitters rated at a maximum output power not listed above, the recommended separation distance D in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

Recommended Separation Distances Between Portable and Mobile RF Communications Equipment and this device. This device and system are NOT Life-Supporting

This device is intended for use in the electromagnetic environment in which radiated disturbances are controlled. The customer or user of this device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF Communications Equipment and the device as recommended below, according to the maximum output power of the communications equipment.

Maximum Output Power (Watts)	Recommended Separation Distances for the device (meters)		
	150 kHz to 80 MHz $D=(1.1667)\sqrt{P}$	80 to 800MHz $D=(1.1667)\sqrt{P}$	800 MHz to 2.5 GHz $D=(2.3333)\sqrt{P}$
0.01	0.11667	0.11667	0.23333
0.1	0.36894	0.36894	0.73785
1	1.1667	1.1667	2.3333
10	3.6894	3.6894	7.3785
100	11.667	11.667	23.333

For transmitters rated at a maximum output power not listed above, the recommended separation distance D in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.





Note 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

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©2011 DeVilbiss Healthcare LLC 4.11
are registered U.S. Patent and Trademark Office
DeVilbiss Healthcare Ltd
DeVilbiss Healthcare Pty. Limited
DeVilbiss Healthcare SAS
DeVilbiss Healthcare GmbH