





## Bravo™ Nasal Pillow System Cleaning and Disinfecting Instructions

The Bravo is intended for use by multiple patients in an institutional setting provided correct disinfection procedures, as listed below, are followed between patients.

These instructions should only be followed by trained healthcare personnel. Please see the Instructions for Use included with the Bravo for further information, including warnings, cautions, and general directions.



### Cleaning and Disinfecting Procedure

1. Remove the pillows from the reservoir. Disassemble the reservoir and the headgear from the Y-circuit.
2. Wash all components with mild soap and warm water. Do not use alcohol or petroleum products or strong cleaners.
3. Make sure that the exhalation ports on the reservoir are free of any obstructions.
4. Rinse thoroughly (3) three times with clear water.
5. If any residue remains on any components after cleaning, discard entire assembly and replace.
6. Prepare a solution of Cidex OPA (*ortho*-Phthalaldehyde solution) per manufacturer's instructions. Verify strength of Cidex OPA solution using Cidex OPA test strips.
7. Submerge all components (except headgear) into the prepared Cidex OPA solution for twelve (12) minutes.
8. Remove components. Shake excess solution from interior surfaces of components.
9. Rinse assembly three (3) times with clear water. Shake excess water from interior surfaces of components. Allow all components to completely air-dry, away from direct sunlight.
10. Carefully inspect each component for any wear or damage. If you see any cracks in the plastic or any tears or hardening of the silicone components, do not use the device. Discoloration is normal after several washings and should not be considered sign of wear. Contact your provider for a replacement.

Technical Support: 800-200-9842, Ordering: 866-216-7337

[WWW.BRAVOMASK.COM](http://WWW.BRAVOMASK.COM)

BRAVO is a trademark of RespCare Inc., Coconut Creek, FL 33073 \*PATENTS PENDING\*