## EN

DeVilbiss® Nonin XPOD Oximeter & Cable Kit  
Instruction Guide, Model DV5M-697

## ES

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**TABLE OF CONTENTS**

Important Safeguards .....	2
Introduction.....	3
Parts and Accessories.....	4
Connecting the Oximeter .....	4
Attaching the Sensor .....	4
Starting Treatment.....	4
Periodic Inspection.....	4
Cleaning the Oximeter.....	4
Specifications .....	5
Declaration of Conformity.....	5

**IMPORTANT SAFEGUARDS**

When using electrical products, basic safety precautions should always be followed. Read all instructions before using. Important information is highlighted by these terms:

- DANGER** Urgent safety information for hazards that will cause serious injury or death.
- WARNING** Important safety information for hazards that may cause serious injury.
- CAUTION** Information for preventing damage to the product.
- NOTE** Information to which you should pay special attention.

**READ ALL INSTRUCTIONS BEFORE USING.  
SAVE THESE INSTRUCTIONS.**

**▲ DANGER**

- **Electric shock hazard - Do not use while bathing.**
- **Electric shock hazard - Do not submerge in water or any other liquid.**
- **Electric shock hazard - Do not attempt to open or remove enclosure; there are no user-serviceable internal components. If service is required, return it to a qualified DeVilbiss provider or authorized service center.**

## WARNING

- **This oximeter is an accessory to the DV5M SmartLink Module. Only connect to this approved product.**
- **Do not use the XPOD oximeter in an MRI environment.**
- **Explosion Hazard– Do not use the XPOD oximeter in an explosive atmosphere or in the presence of flammable anesthetics or gases.**
- **Strangulation Hazard– As with all medical equipment, carefully route cables and connections to reduce the possibility of entanglement or strangulation.**
- **Inspect the sensor application site at least every 6 to 8 hours to ensure correct sensor alignment and skin integrity. Patient sensitivity may vary due to medical status or skin condition. Discontinue use of adhesive tape strips if the patient exhibits an allergic reaction to the adhesive material.**
- **Pacemaker Patients – Rate meters may continue to count the pacemaker rate during occurrences of cardiac arrest or some arrhythmias. Do not rely entirely upon rate meter alarms. Keep pacemaker patients under close surveillance.**

## CAUTION

- To prevent product damage, do not submerge device in either water or any other fluid and do not use caustic or abrasive cleaning agents.
- Loss of monitoring can result if any objects hinder the pulse measurement. Ensure that no blood flow restrictors (e.g., blood pressure cuff) hinder pulse measurements.
- Use only NONIN-manufactured PureLight® pulse oximeter sensors. These sensors are manufactured to meet the accuracy specifications for NONIN pulse oximeters. Using other manufacturers' sensors can result in improper pulse oximeter performance.
- This device does not meet defibrillation-proof requirements per IEC 60601-1: 1990, clause 17.h. To prevent product damage, do not defibrillate a patient while attached to the XPOD oximeter.
- This device should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the device should be observed carefully to verify normal operation.
- The accuracy of the SpO<sub>2</sub> measurement may be affected if the total cable length (including extension cables) is greater than 3 meters.
- The use of accessories, sensors, and cables other than those specified by NONIN (except transducers and cables sold by NONIN as replacement parts for internal components) may result in increased emission and/or decreased immunity of this device.
- This pulse oximetry system is designed to determine the percentage of arterial oxygen saturation of functional hemoglobin. Significant levels of dysfunctional hemoglobin, such as methemoglobin, might affect the accuracy of the measurement.
- The XPOD oximeter has motion tolerant software that minimizes the likelihood of motion artifact being misinterpreted as good pulse quality. In some circumstances, however, the XPOD oximeter may still interpret motion as good pulse quality. This covers all available outputs.
- This equipment complies with International Standard EN 60601-1-2:2001 for electromagnetic compatibility for medical electrical equipment and/or systems. This standard is designed to provide reasonable protection against harmful interference in a typical medical installation. However, because of the proliferation of radio frequency transmitting equipment and other sources of electrical noise in healthcare and other environments, it is possible that high levels of interference due to close proximity or strength of a source might disrupt the device's performance. Medical electrical equipment needs special precautions regarding EMC, and all equipment must be installed and put into service according to the EMC information specified in this document. Portable and mobile RF communications equipment can affect medical electrical equipment.

For more information about required safety and regulatory requirements for medical devices, refer to EN865 and IEC 60601-1. Additional safety information can be found in the labeling provided with each Nonin sensor.

## INTRODUCTION

Your clinician may have supplied an oximeter for use with the DV5M SmartLink system. Connecting an oximeter to the SmartLink allows collection and storage of oximetry data (oxygen saturation and pulse rate).

When using the optional oximeter with the SmartLink, follow the directions below for connecting the oximeter. The only oximeter compatible with this product is the Nonin XPOD Oximeter. Use only the Nonin 7000A, 8000AA, or 8000J sensors. Once the oximeter is plugged into the SmartLink, the oximeter will be powered and ready for use.

## PARTS AND ACCESSORIES

The DV5M-697 Oximeter Kit contains the following items:

1. 9054D-693 Nonin XPOD Oximeter
2. 9054D-695 Nonin Sensor (8000J reusable FlexSensor)

## CONNECTING THE OXIMETER (Figure 1)

1. Plug the oximeter into the bottom circular connector on the side of the DV5M SmartLink Module.
2. Connect the supplied sensor to the oximeter. Match the two arrows on each end of the cables when inserting the oximeter to the sensor.
3. When an oximeter is first connected and the sensor placed onto the finger, a visual notification will be displayed after about ten (10) seconds to signal the presence of new oximetry (SpO<sub>2</sub>) data. SpO<sub>2</sub> will be displayed on the CPAP LCD.
4. Pulse Rate (PR) and oximetry (SpO<sub>2</sub>) data is available on the SmartLink menus.

Connecting the Oximeter to the DV5M SmartLink Module.

A display screen on the LCD display shows real-time SpO<sub>2</sub> and Pulse Rate (updated once every second). Both values will be displayed as “\*\*\*” if the oximeter is not connected or if valid data is not available. This feature is to confirm that the oximeter is connected and working properly. If the “\*\*\*” persists, check connections and verify proper sensor placement and consult the sensor instructions for troubleshooting.

If the sensor is removed from the finger or the signal is lost for more than 10 seconds, the unit will begin to display a visual notification every 8 seconds until the oximeter is re-secured or sensor is replaced onto the finger or a timeout threshold of one minute is reached.

If the oximeter is unplugged from the SmartLink, no notification will occur and the unit will cease logging oximetry data.

The DV5M SmartLink is capable of storing 72 hours of oximetry data (9 days at 8 hours per day).

## ATTACHING THE SENSOR

For instructions on attaching the sensor, refer to the respective pulse oximeter sensor package insert.

**NOTE**—The sensor will not function properly if there is nail polish on your fingernails.

**NOTE**—Check that the sensor is not damaged. If damaged, contact your equipment supplier.

**NOTE**—Do not attach the sensor to the oximeter until you are ready to start using the device.

## STARTING TREATMENT

Start treatment as described in your IntelliPAP™ or SleepCube™ Instruction Guide. The SmartLink will automatically collect and store treatment data.

## PERIODIC INSPECTION

The XPOD oximeter and sensor should be inspected before and after each use. Inspect sensor wires and cables for breaks in the insulation or visible damage. Replace any sensors or cables that have visible damage. Follow the cleaning instructions for the various components of the oximeter system.

## CLEANING THE PULSE OXIMETER

Clean the XPOD oximeter separately from the sensor. For instructions on cleaning pulse oximeter sensor, refer to the respective pulse oximeter sensor package insert(s).

**CAUTION**—To prevent product damage, do not submerge device in either water or any other fluid and do not use caustic or abrasive cleaning agents.

Clean the oximeter components with a soft cloth dampened with isopropyl alcohol. Do not pour or spray any liquids onto the XPOD oximeter and do not allow any liquid to enter any openings in the device. Allow the oximeter components to dry thoroughly before reusing.

## SPECIFICATIONS

Oxygen Saturation Range (SpO <sub>2</sub> ).....	0 to 100%
Pulse Rate Range .....	18 to 300 Pulses Per Minute
Measurement Wavelengths .....	Red - 660 Nanometers @ 3 mW nominal Infrared - 910 Nanometers @ 3 mW nominal

### Accuracy<sup>1</sup>

#### SpO<sub>2</sub> (70-100%) (±1 SD<sup>2</sup>):

No Motion.....	±2 digits
Motion .....	±2 digits
Low Perfusion.....	±2 digits

#### Pulse Rate:

No Motion (18 - 300 BPM).....	±3 digits
Motion (40 - 240 BPM) .....	±5 digits
Low Perfusion (40 - 240 BPM).....	±3 digits

SpO<sub>2</sub> and Pulse Rate display update rate ..... 1 Hz

### Temperature:

Operating.....	32°F to +122°F (0°C to +50°C)
Non-Operating .....	-4°F to +122°F (-20°C to +50°C)

### Humidity:

Operating.....	10 to 90% non-condensing
Non-Operating .....	10 to 95% non-condensing

Power Draw..... 60 mW - typical operation

Dimensions..... 2.1" x 0.8" x 0.6" (53 x 20 x 15mm)

Weight ..... 0.165 lbs. (75g) (including 6' cable and connector)

Sensors ..... Designed to use Nonin sensors only

Equipment classification with respect to protection from electric shock..... Class II

Degree of protection from electric shock..... Type BF Applied Part

Warranty..... One-year limited

## DECLARATION OF CONFORMITY

**Manufacturer:** DeVilbiss Healthcare  
**Address:** 100 DeVilbiss Drive, Somerset, PA 15501  
**Product Designation:** CPAP Accessory  
**Type, Model:** DV5M SmartLink Module  
 For DeVilbiss IntelliPAP or SleepCube CPAP

We herewith declare that the above mentioned product complies with the requirements of EC Directive 93/42/EEC and the following:

**Class:** IIa, Rule 2

### Quality System

**Standards Applied:** ISO 13485:2003

**Notified Body:** TUV NORD

**MDD:** Annex II Applied

**Safety Standards Applied:** IEC 60601-1:1988 + A1 1995

IEC 60601-1-4:2000

ISO 14971:2000

**EMC Compliance To:** IEC 60601-1-2:2001

**European Contact:** Sunrise Medical Ltd.  
 High Street  
 Wollaston West Midlands DY8 4PS  
 England  
 44-138-444-6688

<sup>1</sup> All accuracy specifications are results determined by induced Hypoxia studies on healthy adult volunteers using the 8000AA Finger Clip Sensor.

<sup>2</sup> Standard Deviation (SD) is a statistical measure: up to 32% of the readings may fall outside these limits.

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A display screen on the LCD display shows real-time SpO<sub>2</sub> and Pulse Rate (updated once every second). Both values will be displayed as “\*\*\*\*” if the oximeter is not connected or if valid data is not available. This feature is to confirm that the oximeter is connected and working properly. If the “\*\*\*\*” persists, check connections and verify proper sensor placement and consult the sensor instructions for troubleshooting.

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**Sunrise Medical Pty. Limited**
Nonin Medical, Inc.