DEVILBISS®

DeVilbiss® PulseDose® LT Compact Conserving Device

Model PD1000 Instruction Guide

CAUTION— Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

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FRANÇAIS 40

TABLE OF CONTENTS

Important Safeguards	5
Introduction	8
Use Times	9
Important Parts of Your PulseDose Compact Conserving Device	10
Operating Instructions	12
Typical Questions and Answers	14
Care & Maintenance	16
Specifications	17
Troubleshooting	18
Important Information	20

IMPORTANT SAFEGUARDS

The information contained in this guide is intended to assist in the safe operation of the equipment and to ensure maximum benefit is achieved.

This product is to be used only to deliver medical grade (U.S.P.) oxygen and only with a physician's prescription.

READ ALL INSTRUCTIONS BEFORE USING.

SAVE THESE INSTRUCTIONS.

Oxygen supplied by this equipment is not to be considered life-supporting and must not supply anything other than medical grade (U.S.P.) oxygen.

When using electrical products, basic safety precautions should always be followed. Read all instructions before using.

Important information is highlighted by these terms:

DANGER– Urgent safety information for hazards that will cause serious injury or death.

WARNING– Important safety information for hazards that might cause serious injury.

CAUTION– Information for preventing damage to the product.

NOTE– Information to which you should pay special attention.

 **ATTENTION**– Consult accompanying documents.

 **DANGER**– **NO SMOKING**

 Type BF applied part

⚠ DANGER

To reduce the risk of fire, burns, or injury to persons:

Oxygen, though non-flammable, vigorously supports and accelerates burning of any flammable material. If you know or suspect oxygen has escaped other than through normal operation, open doors and windows to ventilate the area.

1. DO NOT SMOKE WHILE USING YOUR DEVILBISS OXYGEN EQUIPMENT.
Keep matches, cigarettes, burning tobacco, or candles away from the area where the system is being stored or operated.
2. Avoid creation of any spark near oxygen equipment. This includes sparks from static electricity created by any type of friction.
3. Keep the equipment at least seven feet away from radios, television sets, window air conditioners, fans, electric razors, hair dryers, and all other electrical appliances.
4. Keep the equipment away from heat sources, electric or gas heaters of any kind, fireplaces, or stoves.
5. Keep all flammable materials or petroleum-based products away from the equipment.
6. Never attempt to lubricate the equipment.
7. Never use aerosol sprays near the equipment.

To prevent high concentrations of oxygen:

1. Keep the equipment in a well-ventilated area.
2. Do not carry equipment under a coat or any form of clothing.
3. Turn off oxygen supply by closing the cylinder valve when not in use.

⚠ WARNING


To reduce the risk of injury:

1. Keep all units away from children. Do not allow unauthorized or untrained individuals to operate the equipment. Never tamper with or try to repair the equipment yourself. If you have any questions or suspect your equipment is not operating properly, contact your oxygen provider.
2. Do not immerse in liquids or subject device to harsh conditions.
3. Do not use in temperatures greater than 104°F (40°C) or below 41°F (5°C).
4. Do not use with other equipment (i.e. humidifier, nebulizer, etc.) when in PulseDose delivery mode.

PHYSICIAN'S NOTES

1. Do not use with patients who breathe below 6 Breaths Per Minute (BPM) or above 40 BPM.
2. Do not use with patients who consistently fail to trigger equipment (i.e. mouth breathing with closed soft pallets).
3. Verify patient is getting adequate PaO₂ or SaO₂ levels in PulseDose delivery.
4. Use only standard nasal cannula with PulseDose delivery. Do not use pediatric (low-flow) nasal cannula or mask with PulseDose delivery.
5. A mask or any nasal cannula can be used with continuous flow delivery.
6. The standard device contains a preset 2 LPM continuous flow backup.
The service manual describes how to change the cannula fitting to obtain a 3, 4, 5 or 6 LPM continuous flow backup.

NOTE– This conserving device is designed to fit onto 870-style or “pin-index” cylinders as pictured.

INDICATIONS FOR USE

The DeVilbiss PD1000 Compact Conserving Device is intended as a delivery device for medical-grade oxygen from high-pressure oxygen cylinders. This is an ambulatory device, which allows patients to ambulate longer than they would with a continuous flow regulator on the same cylinder.

HOW PULSEDSE WORKS

PulseDose dramatically extends the use time from a supply of oxygen, offering increased mobility with improved comfort and increased efficiency. The reliability and safety of PulseDose oxygen delivery has been proven effective in clinical testing as well as through independent tests performed by physicians and respiratory therapists.

What is PulseDose? The concept is based on the fact that the normal breathing pattern is inhalation for 1/3 of the time, and exhalation about 2/3 of the time. At 20 BPM the oxygen delivered in continuous flow, assuming inspiration is 1/3 of the breathing cycle, would be 16.5 cc/LPM. As a result, PulseDose extends the use time of an oxygen system by an average of 3:1. PulseDose senses the start of inhalation and instantly releases a short “pulsed” dose at the very beginning of the breathing cycle. Since all of the “pulsed” oxygen finds its way deep into the lungs, less oxygen is required to accomplish the same effect than with traditional continuous flow oxygen systems. This means that a PulseDose oxygen system will last two to four times longer than a continuous flow oxygen system, yet still provide the same therapeutic benefit.

Because oxygen is released only during inhalation, the constant flow of oxygen into the nostrils is eliminated. Many users find PulseDose oxygen delivery more comfortable than continuous flow delivery systems. The short “pulse” of oxygen delivered during inhalation is almost undetectable, and the humidity in the room air helps maintain a normal level of moisture in the nasal cavity. This greatly reduces the discomfort of dehydration associated with a continuous flow oxygen system.

Because PulseDose responds to each individual’s breathing patterns, the use time will vary for each individual depending on the PulseDose prescription rate and the breath rate. The following chart shows the theoretical ambulatory ranges for DeVilbiss PulseDose products.

NOTE– All ambulatory ranges are calculated assuming a breath rate of 20 breaths per minute in PulseDose (PD) mode.


USE TIMES SHOWN IN HOURS									
Delivered Volume cc’s:	16.5	24.75	33	41.25	49.5	66	82.5	99	
Flow Rate:	1	1.5	2	2.5	3	4	5	6	Mode
Walkabout Mini M-6 Cylinder 164 Gaseous Liters	2.7	1.8	1.4	1.1	.9	.7	.6	.4	CF
	8.3	5.5	4.1	3.3	2.8	2.1	1.7	1.4	PD
Hideaway ML-6 Cylinder 170 Gaseous Liters	2.8	1.9	1.4	1.1	.9	.7	.6	.5	CF
	8.6	5.7	4.3	3.4	2.9	2.1	1.7	1.4	PD
Walkabout 1 “C” Cylinder 240 Gaseous Liters	4.0	2.7	2.0	1.6	1.3	1.0	.8	.7	CF
	12.1	8.1	6.1	4.9	4.0	3.0	2.4	2.0	PD
Walkabout 2 “D” Cylinder 415 Gaseous Liters	6.9	4.6	3.5	2.8	2.3	1.7	1.4	1.2	CF
	21.0	14.0	10.5	8.4	7.0	5.2	4.2	3.5	PD
Walkabout 3 “E” Cylinder 682 Gaseous Liters	11.4	7.6	5.7	4.6	3.8	2.8	2.3	1.9	CF
	34.4	23.0	17.2	13.8	11.5	8.6	6.9	5.8	PD

CF=Continuous Flow

PD=PulseDose

Specifications subject to change without notice. This chart is intended to be used only as a guide.

Cylinders vary in gaseous liter capacity by manufacturer which may result in varying use times.

1. **Knob**– This knob is used to attach the conserving device to the cylinder.
 2. **PulseDose Indicator**– Either green or red light illuminates each time the unit pulses oxygen.
 - **Normal Battery Indicator**– A flashing green light indicates that there is sufficient battery power.
 3. **Low Battery Indicator**– A flashing red light indicates that there are 4 - 8 hours of battery life remaining. The low battery indication time may be reduced when using NiMH batteries.
 - **Change Battery Indicator**– A constant red light indicates that the battery should be changed immediately. Unit can be used only in continuous flow mode until new batteries are installed.
 4. **Cannula Fitting**– Use this fitting to attach the cannula to your PulseDose conserving device.
 5. **Rotary Selector**– When this rotary switch is set to “OFF,” the unit is not using battery power and will not pulse. When this switch is set to one of the numbers, the unit is on and awaiting inspiration through the nasal cannula at which time it will dose on every breath. The volume of the oxygen delivered varies according to which prescription flow setting is chosen. The final setting on the rotary switch is “CF”; this is the continuous flow position. In this position oxygen will flow from the cannula fitting at the default continuous flow rate.
 6. **Oxygen Contents Gauge**– Indicates the remaining pressure in the oxygen cylinder. When this gauge falls into the red section, you should switch to a new cylinder.
 7. **Battery Door**– Use only a standard “AA” Alkaline or NiMH batteries.
- NOTE**– Continuous flow mode is not powered by the batteries and can be used regardless of the battery level. In the event of a device failure or dead batteries, the user must manually switch the unit to continuous flow for delivery of oxygen. The device will not automatically switch to continuous flow. The oxygen cylinder will not last as long in continuous flow mode as it would in PulseDose mode. Unless there is a problem with the unit, such as dead batteries, the unit should be used in PulseDose mode.

INSERTING A BATTERY INTO THE PULSEDSE COMPACT CONSERVING DEVICE

1. Open the battery door.
2. Insert 2 "AA" Alkaline or NiMH batteries.
3. Close the battery door.

NOTE– When changing batteries, first turn the rotary selector to the "OFF" position.

ATTACHING YOUR PULSEDSE COMPACT CONSERVING DEVICE TO THE CYLINDER

1. Loosen the knob.
2. Carefully lower the conserving device over the post of the cylinder.

CAUTION– Alignment pins can damage sealing surfaces of the post increasing the chance of leakage.

3. Align the pins in the conserving device to the holes in the cylinder post as you would a standard regulator.
4. Hand-tighten the knob until the conserving device is secure.

USING YOUR PULSEDSE COMPACT CONSERVING DEVICE

1. Open the cylinder.
2. Attach the standard nasal cannula to the conserving device and to your nose and face. Oxygen tubing up to 35 feet long may be used in PulseDose delivery mode.
3. Turn the rotary selector to the prescribed flow setting.
4. Breathe normally, the conserving device will deliver a bolus of oxygen at the leading edge of inspiration on every breath up to 40 breaths per minute.
5. When you are finished using the PulseDose conserving device, turn your cylinder to the closed position, and rotate the rotary selector to the "OFF" position.

⚠WARNING

To prevent injury from cylinders tipping over, do not use cannula tubing lengths over 10 feet with small compressed oxygen cylinders. Unattended cylinders should be secured in a cylinder stand.

NOTE– A mask should not be used in the PulseDose delivery mode as it may not fit to the face well enough to allow the conserving device to sense inhalation efforts. Also, the therapeutic effect of PulseDose would not be realized, as the dose of oxygen would be diluted in the mask prior to inhalation.

NOTE– A pediatric or low-flow cannula should not be used in PulseDose delivery mode. The reduced diameter of the cannula causes too much back pressure and will affect the oxygen volume delivered.

NOTE– PulseDose delivers oxygen in a very short "puff." It does not deliver oxygen during the entire inhalation. The length of time that PulseDose delivers oxygen will not vary from breath to breath. The time is set in correlation to the oxygen dosage set on the conserving device (patient's prescription setting).

NOTE– PulseDose is designed to prevent the delivery of pulses more than every 1 1/2 seconds. If the breath rate is greater than 40 BPM, this feature prevents delivery of excessive oxygen by not dosing on every breath.

NOTE– If using NiMH batteries, carefully monitor when the low battery indicator flashes red. It is recommended that a spare fully-charged set of alkaline batteries be kept in reserve and installed when the change battery indicator is constant red.

Q. How does PulseDose work? How does it know when I'm inhaling?

A. When inhaling, your diaphragm moves down and causes a drop in pressure in the lungs. Air flows in through the nose and mouth to equalize the pressure. This negative pressure is also present at the nose and mouth during inhalation. This pressure signal travels through the nasal cannula to a pressure sensor in the PulseDose conserving device. An electronic circuit then opens an electrical valve to deliver a precisely metered dose of oxygen. When the valve is closed, the sensor is ready to detect the next inhalation.

Q. The pulse seems so short. Am I really getting enough oxygen?

A. Yes. PulseDose delivers a precise burst of oxygen at a relatively high flow rate at the leading edge of each inhalation. This assures that the oxygen delivered flows deep into the lungs for maximum benefit. PulseDose requires less oxygen to deliver the same therapeutic benefit as continuous flow oxygen delivery.

Q. I can't hear the pulse. Is PulseDose working?

A. If the pulse can't be heard, simply look at the green PulseDose indicator to see that the device is being triggered by inhalation. For further assurance, hold the end of the cannula in front of your lips while inhaling through your mouth and feel the pulse. PulseDose does not monitor the supply of oxygen. Remember to check the oxygen contents gauge periodically to verify that there is an adequate oxygen supply. If the oxygen supply runs out, the green PulseDose indicator light will continue to illuminate, indicating that the conserving device is being triggered by inhalation.

Q. Why can't I use a cannula which is longer than 35 feet?

A. The PulseDose triggering is not significantly affected by the cannula length, but the delivery of oxygen is affected. If the cannula is longer than 35 feet, the pulse of oxygen is delayed. Remember the therapeutic moment during the inhalation cycle. If the oxygen is not delivered during this time, the benefits will not be realized.

Q. I've always used humidifiers with oxygen. Should I use a humidifier with PulseDose?

A. No. PulseDose is not able to sense inhalation through the water in the humidifier. Also, many patients find that humidification is not necessary with PulseDose. They find that PulseDose improves comfort because it delivers a very small amount of oxygen during the early part of inhalation, while the rest of the inhalation is composed of normal room air.

Q. When I'm breathing faster, I don't get a pulse with each breath. Don't I need a dose every time?

A. Because PulseDose breathes with the patient, it has an upper limit (40 Breaths Per Minute) that keeps you from getting too much oxygen. When breathing slowly, you receive a dose with every breath. As breath rate increases (up to 40 BPM) PulseDose still delivers a dose with every breath. At this point, you are getting more oxygen per minute because each pulse delivers the same amount of oxygen with each breath while the number of breaths has increased. With continuous flow oxygen, the oxygen delivered is constant. As you breathe faster, the enrichment of inhalations actually decreases because each breath is being diluted with a greater amount of room air.

CARE & MAINTENANCE

The DeVilbiss PulseDose compact conserving device should be kept clean and free from moisture and dust. The device should be protected from extreme temperatures. Clean the device periodically by wiping it with a dry, lint-free cloth. Do not clean the unit with a solvent based cleaning solution. Avoid dropping the conserving device or placing it in a position where it could topple or fall since this can damage the device. Whenever possible, use a padded carrying bag such as the Hideaway or Walkabout bags to carry the conserving device and cylinder. This will help to protect the conserving device in case of a fall. Avoid getting fluids or debris such as sand or dirt inside the device.

CAUTION– To prevent product damage, do not expose the PulseDose conserving device to water.

SPECIFICATIONS

Weight	14.7 ounces 16.3 ounces with battery
Dimensions	4.75”L x 3.4”W x 2.8”H (12.06 cm L x 8.64 cm W x 7.11 cm H)
Power Supply	(2) Standard “AA” alkaline or NiMH batteries.
Power Requirements.....	Batteries other than alkaline or NiMH are not recommended due to the capacity needed for operation and battery life of the unit.
Operating Temperature Range.....	5° to 40°C (41° to 104°F)
Operating Pressure Range	500 to 2250 PSIG (34 to 155 bar) tank pressure
Operating Atmospheric Conditions.....	500 to 1020 millibar
Operating Humidity Range	0 to 95% R.H., non-condensing
Storage and Transportation Temperature Range.....	-20° to 60°C (-4° to 140°F)
Storage and Transportation Humidity Range	Up to 95% R.H., non-condensing
Degree of Protection Against Ingress of Liquids	NONE
Degree of Protection Against Electric Shock.....	TYPE BF applied part
Modes of Operation.....	Continuous / Pulsed
Approval Body And Safety Standard	IEC 601-1; CAN/CSA-C22.2 No 601.1-M90
US Patents	4,519,387; 5,755,224; 4,457,303

⚠ WARNING

Do not attempt to open the device for maintenance or repair. The conserving device contains no user-serviceable parts. Contact your oxygen provider if service is required.

Symptoms	Possible Causes	Remedies
Oxygen is not being delivered even though the PulseDose indicator is flashing every time I breathe.	<ol style="list-style-type: none"> Oxygen supply is empty. Oxygen supply is not turned on. 	<ol style="list-style-type: none"> Check contents indicator on the device. If empty, switch cylinders. Open the compressed oxygen cylinder valve by following the directions given by your service representative.
Use times are different from those stated in the literature.	<ol style="list-style-type: none"> PulseDose responds to your breath rate. Your breath rate may vary, which causes the operation time to vary. Leak in system 	<ol style="list-style-type: none"> PulseDose is operating correctly. Check connection to cylinder. It may need a new regulator seal.

Symptoms	Possible Causes	Remedies
PulseDose will not pulse	<ol style="list-style-type: none"> Cannula is not adjusted properly. Unit is not turned on. Batteries discharged or not installed. Mouth breathing with closed soft palates. Unit did not reset while changing batteries (red light stays on). 	<ol style="list-style-type: none"> Check all cannula connections to make sure they are tight, and adjust the cannula to fit comfortably in your nose. Ensure tubing is not kinked. Turn the rotary selector to the appropriate setting. Install new batteries. Breathe through the nose (cannula). Turn the unit off and back on using the rotary selector.
PulseDose works fine for a couple of minutes, then sensitivity seems to drift and may stop working altogether.	Using pediatric cannula or any cannula that restricts continuous flow capacity of 10 lpm.	Replace with standard nasal cannula.

PHYSICIAN INFORMATION

Name: _____

Address: _____

Telephone: _____

Emergency Telephone: _____

PRESCRIPTION INFORMATION

Patient's Name: _____

Flow Setting (LPM): _____

SET-UP INFORMATION

Name of Person Setting Up: _____

OXYGEN PROVIDER

Emergency Telephone Number: _____

This instruction guide was reviewed with me and I have been instructed on the safe use and care of the DeVilbiss PulseDose oxygen conserving device.

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©2002 Sunrise Medical 5.02
Sunrise Medical Ltd.
Sunrise Medical Pty. Limited