

## Nebulizers, Masks and Accessories

### Efficient solutions for delivering aerosol therapy

DeVilbiss introduces to its compressor accessory line several new products designed to **provide efficient and affordable solutions** to common issues. For pediatric patients, the new Pediatric Puppy Character Mask makes treatment time fun, and its **customizable fit** ensures that children remain comfortable throughout the treatment. This child-friendly mask **features a swivel connector for aerosol treatment in almost every position**, while its custom vent system **keeps aerosol out of the child's eyes**. For infants, the patented PediNeb™ system provides nasal delivery of aerosol while the pacifier soothes the patient.



Pediatric Puppy  
Character Mask



VixOne  
PediNeb



VixOne Disposable  
Nebulizer



0044  
These products are latex-free

The VixOne nebulizer, available in reusable and disposable models, is the **economical, high-performance choice in small-volume nebulizers**. Its ideal particle size production allows for high-volume lung deposition, resulting in maximum therapeutic effect. The corrugated tubing, standard with every unit, provides a superior delivery system resulting in **18% more medication being delivered to patients than with a traditional open system**.

### DeVilbiss Breathies Bunch™

This unique program targets the youngest patients and is packed with age appropriate information, including a guide to childhood asthma, entertaining yet educational videos, log book, colorful stickers and more.

## VixOne™ Nebulizers

- Produces an aerosol with an MMAD of 2.7 microns
- Respirable Fraction – 84% below 5 microns
- **At least 18% more medication delivered** than with a traditional open system
- Low residual volume (0.7cc) reduces medication waste
- High-capacity medication bowl (>10cc) accommodates varying therapy protocols
- Nebulizes in upright position or at angles up to 45°, allowing patients to **receive treatment in their most comfortable position**
- **Anti-spill tee design** provides secure seal with simple ¼ turn assembly design

## Aerosol Masks



Pedi Aerosol Mask



Pediatric Puppy Character Mask



Adult Aerosol Mask

- Soft, pliable, latex-free masks conform comfortably around a wide array of facial contours.

## PediNeb™ Pacifier



The PediNeb was specially designed to enhance the delivery of medication to infants and toddlers while reducing the apprehension experienced by patients.

- Pacifier calms the child during aerosol treatment
- Medication port directs medication to the nose and reduces medication waste associated with the “blow-by” technique
- Easy to operate, clean and maintain
- Patented design was developed by a pediatric nurse practitioner

## Accessories

ITEM	DEVILBISS PART #	QUANTITY
MicroMist Disposable Nebulizer	4650D-621	50 per case
VixOne Disposable Nebulizer	3655D-621	50 per case
VixOne Reusable Nebulizer	900D-610	10 per case
VixOne Disposable with PediNeb®	3655PN-612	12 per case
Pediatric Puppy Character Mask	DL1050	50 per case
Adult Aerosol Mask	AAM12	12 per case
	AAM20	50 per case
Pediatric Aerosol Mask	PAM10	50 per case
Peak Flow Meter	600-112	12 per case
Breathies Bunch Compliance Kit	DBB-KIDS	Each



DeVilbiss Healthcare | 100 DeVilbiss Drive | Somerset, PA 15501 USA  
800.338.1988 | 814.443.4881 | www.DeVilbissHealthcare.com