

Bravo[®] Nasal Pillow System

A Brilliant Design for Nasal CPAP Interfaces

The Bravo is a **unique development in nasal pillow interfaces**, designed to promote patient comfort, allow for a clear field of vision, and **provide the patient multiple wearing options**. With its close fitting reservoir and secure headgear, this **lightweight interface** can be worn over the head or on the chest, and will not interfere with wearing glasses.

The Bravo features the versatility of **interchangeable nasal pillows** and Adaptive-Flow™ technology resulting in an **extremely quiet sleep experience**. The unprecedented Active-Swivel™ feature allows for multiple wearing positions resulting in the **highest level of patient comfort**.



Bravo
Nasal Pillow System
(BRV600)

Performance Features

- Adaptive-Flow technology facilitates **ultra-quiet operation at high pressures**
- **Multiple wearing options** to match the patient's sleep positioning
- Patient can wear glasses while wearing the interface
- **Limited points of contact on face** and no air leaks near eyes
- No forehead or overhead contact

Design Features

- 100% latex free
- **One package contains all of the components** required to properly size the patient, eliminating the hassle and waste of opening multiple products to properly fit the patient
- **Interchangeable nasal pillow inserts** can be mixed, matched and cut to fit each nostril independently
- Nasal pillows designed to insert into the nostrils at an angle for a more natural and comfortable fit
- Ergonomic reservoir design contoured to best fit facial profile
- **Four-point adjustable headgear** provides an optimal seal and increased support
- Sturdy headgear can be applied to head and removed without undoing any straps which provides a consistent fit
- Patient Instructional DVD and carrying case included with each unit
- To view the video please visit: www.BravoMask.com

Reimbursement information

Reimbursement Information

Product Description	Code Description	HCPCS Coding	Frequency	Allowable Range
BRAVO PILLOWS	Pillows for nasal application device	A7033	Two pairs per month	\$24 -\$28
BRAVO APPLICATION	Nasal interface (mask or cannula type) used with positive airway pressure device, with or without head strap	A7034	One per 3 months	\$99 -\$117
HEADGEAR	Headgear used with positive airway pressure device	A7036	One per 6 months	\$33 -\$39

Accessories and Replacement Parts

PART NUMBER

PRODUCT DESCRIPTION

BRV600	Complete Bravo System (Includes: small, medium, large Pillows and Headgear)
BRV630	Bravo Headgear
BRV625	Bravo Pillows Large
BRV623	Bravo Pillows Medium
BRV621	Bravo Pillows Small
X2009	Chin Strap



DeVilbiss Healthcare LLC | 100 DeVilbiss Drive | Somerset, PA 15501 USA
800.338.1988 | 814.443.4881 | www.DeVilbissHealthcare.com