

DeVilbiss FlexSet® Nasal CPAP Mask

One Mask, So Many Options

With so many nasal CPAP mask options on the market, it is cost prohibitive to stock all of them. The DeVilbiss FlexSet mask provides four different cushion options, allowing you to minimize the masks in stock and **build the mask to suit your patient.**

The FlexSet is amazingly intuitive to use and care for. This innovative interface features a unique forehead adjustment and pivoting design that help to ensure a **secure, leak-free fit on most facial profiles.** Its lightweight design and minimal headgear make the mask **extremely comfortable** for the patient to wear. Additionally, the tubing can be worn straight down or over-the-head for added flexibility.



DeVilbiss FlexSet
Nasal CPAP Mask, Standard
(9354D)

Performance Features

- Patented **adjustable exhaust port** allows patient to direct escaping air away from their chest or bed partner
- One-piece, **self-anchoring headgear** is easy to use
- Dual-pivoting forehead rest and mask bridge for greater comfort and improved fit
- Two ultra-soft Breath-O-Prene® forehead pads included for optimal patient comfort
- Hose connection swivel elbow **rotates 360 degrees** or can be locked in an upright position to run tubing over the head

Design Features

- 100% latex free
- Soft, replaceable cushion provides a **comfortable low pressure seal**
Available in standard and shallow sizes and silicone or gel material
- **Quick disconnect for ease of use** – tubing can be disconnected from the mask swivel
- Minimal headgear offers a secure fit and features a quick release buckle.
- Supplemental oxygen port.
- The FlexSet mask provides a modular and easy-to-use nasal CPAP interface that allows users to build a mask to suit their needs.
- There are currently 4 different cushion options available (see accessory and replacement parts listing below) and each FlexSet mask comes complete with 2 different forehead pad options.
- The FlexSet Mask is dishwasher safe (with cushion, headgear and forehead pad removed).

Reimbursement information

Product Description	Code Description	HCPCS Code	Frequency	Allowable Range
FlexSet Mask	Nasal Application Device (mask or cannula type), used with positive airway pressure device, with or without head strap	A7034	1 per 3 months	\$99.00 - \$117.64 initial set-up
Replacement Cushion	Replacement cushion for Nasal Application Device, each	A7032	2 per 1 month	\$34.45 - \$40.53
Replacement Headgear	Headgear	A7035	1 per 6 months	\$33.79 - \$39.75
Chin Strap	Chinstrap	A7036	1 per 6 months	\$15.47 - \$18.20

Accessory and Replacement Parts

ITEM

FlexSet Mask Standard
 FlexSet Mask Shallow
 FlexSet Mask Gel
 FlexSet Mask Gel Shallow
 FlexSet Silicone Standard Cushion
 FlexSet Gel Standard Cushion
 FlexSet Silicone Shallow Cushion
 FlexSet Gel Shallow Cushion
 FlexSet Headgear
 FlexSet Chin Strap
 FlexSet Pad/Vent

DEVILBISS PART

9354D
 9354S
 9354G
 9354GS
 9352D-601
 9352G-601
 9352S-601
 9352GS-601
 9353D-619
 7351D-647
 9353D-604



Standard and Shallow Gel Cushions
 9352G-601, 9352GS-601



DeVilbiss Healthcare LLC | 100 DeVilbiss Drive | Somerset, PA 15501 USA
 800.338.1988 | 814.443.4881 | www.DeVilbissHealthcare.com