

## HomeCare Suction Unit

### A Fusion of Quality and Innovation

For more than twenty years, DeVilbiss has been delivering the **best portable suction options**. The DeVilbiss 7305 Series Homecare Suction Unit **meets suction and capacity requirements in a small, lightweight design**.

The DeVilbiss 7305 series Homecare Suction Unit is a compact medical suctioning device which has been **designed for reliable, portable operation**. Because of the **small size, light weight and DC operation**, the DeVilbiss 7305 Suction Unit is ideal for providing suction in the home, or on the go with the **DC cord**. The 7305P series with **internal rechargeable battery** also available. **Two collection bottle options** give the choice between the standard **800 ml disposable canister** and optional long-term **reusable 1,200 ml canister**.



HomeCare Suction Unit  
Model 7305P series with  
internal rechargeable  
battery also available

### Performance Features

- High Performance - Vacuum adjustments allow for 50-550 mm Hg and a free flow of 27 liters/minute.
- Adjustable Flow Regulator - Easy-to-turn knob conveniently placed directly below gauge, for easy use and control. Knob has safety lock feature which will not allow the knob to be twisted off.
- The DeVilbiss® Homecare Suction Unit 7305P series meets the RTCA/DO160D Aircraft Standard as well as the ISO-10079.1:1991 Suction Standard
- Increased Battery Capabilities – For use up to one hour. This provides the patient with added security when it is needed most.

## Design Features

- Easy-to-Read Gauge - The unit has conveniently placed gauge for easy reading. This is especially important when precise measurements of vacuum are needed for the pediatric patient.
- The 7305 series is equipped with a switch mode power supply allowing operation on any AC voltage (100-240 VAC, 50/60 Hz).
- Two-year warranty on Suction Unit.
- Compact Design - All of the high-performance capabilities are conveniently placed in a small compact package weighing less than 6.3 lbs. The compact design with the high-performance capabilities makes the 7305 series suction unit the smallest and strongest portable suction unit of its kind.
- Low battery and charge indicator lights on portable unit (7305P-D).

## Accessories and Replacement Parts

### Available Models

**7305D-D Suction Unit, 800cc disposable bottle, tubing, elbow, and AC & DC power cords**

**7305P-D Suction Unit, 800cc disposable bottle, tubing, elbow, AC & DC power cords, and rechargeable battery**

PART NUMBER	PRODUCT DESCRIPTION	PART NUMBER	PRODUCT DESCRIPTION
6305D-611	6' Patient Tubing	7305P-610	Circuit board (7305P series & converted battery units)
7304D-619	12V DC power cord (1ea)	7305D-610	Circuit board (7305D series)
7305D-601	1200 ml Reusable collection bottle (16 ea)	7305P-620	High capacity battery conversion kit (battery, circuit board & battery tray)
7305D-602	800 ml Disposable collection bottle (48 ea)	7304D-636	Vacuum gauge test fixture
7305P-613	Switch mode AC to DC adapter/charger	099HD-614	Hospital grade line cord (120 VAC)
7305P-614	12V DC high capacity rechargeable battery (1 ea)	6710D-609	Line cord for US
7305D-606	Carrying case	7305P-631	Line cord for Continental Europe
7305D-607	Replacement gauge	7305D-603	Filter, 800ml collection bottle, elbow & tubing package
7305D-608	Bacteria filter (non-sterile) (12 pack)	7305D-604	1200ml bottle package

## Technical Specifications

<b>Dimensions</b>	9" x 7.0" x 8.0" (22.9cm x 17.8cm x 20.3cm)
<b>7305D Series</b>	3.8 lbs (1.7 kg)
<b>7305P Series</b>	6.3 lbs (2.9 kg)
<b>Electrical Requirements</b>	100-240 V AC 50/60 Hz .75 A max; 12 V DC, 33 W max
<b>7305P Series Internal Rechargeable Battery</b>	Factory Equipped
<b>7305D Series Internal Rechargeable Battery</b>	Not Factory Equipped (provider-installed option)
<b>7305P &amp; 7305D Series Vacuum Range</b>	80 to 550 mm Hg
<b>Air Flow @ pump inlet:</b>	27 LPM (free flow) typical (may be less when running from internal battery)
<b>7305 D&amp;P Series Collection Bottle Capacity</b>	800 ml (cc) Disposable
<b>7305 D&amp;P Series Collection Bottle Capacity-Optional</b>	1,200 ml (cc) Reusable
<b>7305P Series Warranty</b>	Two-years limited, excluding internal battery and collection bottle
<b>7305D Series Warranty</b>	Two-years limited, excluding collection bottle
<b>Internal Battery Warranty</b>	90-day
<b>7305D &amp; 7305P Series Approvals</b>	IEC 601-1; CAN/CSA-C22.2 No. 601.1-M90; UL 60601-1, EN 60601-1-2
<b>7305P meets RTCA/DO &amp; 160D Approvals for battery operation only; Airline use.</b>	DO-160D - section 21 Category M, ISO-10079.1:1991



DeVilbiss Healthcare | 100 DeVilbiss Drive | Somerset, PA 15501 USA  
800.338.1988 | 814.443.4881 | [www.DeVilbissHealthcare.com](http://www.DeVilbissHealthcare.com)