



DeVilbiss IntelliPAP® Standard CPAP System

You talked. We listened.
Everyone gets what they need.



DeVilbiss IntelliPAP Standard
CPAP System (Model DV51D-HH)

The **DeVilbiss IntelliPAP Standard CPAP System** was **designed with patients and providers in mind**, incorporating many of the recommendations solicited through research to **optimize patient adherence**. IntelliPAP's combination of **comfort, education and adherence tracking with SmartCode®** helps to ensure patient compliance – all in a **highly efficient platform** that makes great business sense to you.

From its **reliable electronics and water resistant housing** to its durable Lexan® humidifier chamber and stainless steel heat transfer plate, IntelliPAP is **built for years of dependable performance and patient satisfaction**. It's just what you would expect from DeVilbiss – the widely recognized leader in respiratory technology.

Performance Features

- **Among the quietest CPAPs** on the market – ONLY 26 dBA!*
- Optional integrated heated humidification system
- **Automatic Altitude Adjustment** up to 9,000 feet
- Exclusive humidifier chamber design **maximizes humidification and minimizes leaks**
- Exceptionally small footprint – **smallest integrated unit** among top-selling CPAP devices
- Delay ramp pressure eases unit to prescribed pressure
- Best-in-industry warranty

Design Features

- Clean, sophisticated design
- Whisper-quiet operation – 26 dBA!*
- Exceptionally small footprint fits easily into the home environment and carry-on luggage
- Includes discreet carry bag that can be easily expanded to hold additional accessories
- Optional integrated heated humidification system easily docks into the bottom of the unit without requiring a separate hose or power cord
- Humidifier reservoir is easy to remove and replace, and includes built-in overflow protection
- Water ingress protection prevents humidifier water from damaging the CPAP – great for travelers
- Patient air tubing extends from back of unit to reduce clutter and enhance visual appeal
- Accommodates various tubing lengths (6 and 10 feet)
- Large, slip-resistant pads on bottom of unit to minimize movement on night table during the night
- Eight-foot power cord length provides greater positioning flexibility – no bulky power brick
- Easy and intuitive set up – three, large buttons control nearly all functions for simplified patient operation
- Large, backlit LCD and easy-to-follow patient menu
- Auto on/off** automatically starts the flow of air after breathing into the mask and stops the flow of air when mask is removed
- Delay pressure ramp eases the machine to the prescribed pressure – allowing the user to fall asleep at a lower pressure (0-45 minutes in 5-minute increments and .5 increments for pressure)
- SmartCode compliance tracking encodes usage data that can be communicated to a clinician and then translated
- Reminders alert user to replace or check parts
- Extra air inlet and fine particle filters included
- Industry-leading 3-year warranty

Patient Education

- User video includes an introduction to CPAP therapy with tips for acclimating to and optimizing therapy, as well as how to care for and use the IntelliPAP system
- Patient instruction guide demonstrates how easy it is to use IntelliPAP, while providing a comprehensive overview of sleep apnea, its treatment and tips for acclimating to therapy.

DeVilbiss IntelliPAP Standard CPAP System includes attractive travel bag (DV51D)



Technical Specifications

Product Dimensions	6.4" H x 6.5" W x 8.4" D	Automatic Leak Compensation	Yes
Product Weight	4.45 lbs with humidifier and PAP	Auto On/Off	Yes (can be enabled or disabled)
Pressure Range	3-20 cm H ₂ O (.5 increments)	Motor Hour Meter	Yes
Universal Power Supply	Yes	Therapy Compliance Meter	Yes
AC Cord Length	8 feet	Warranty	3-year warranty on IntelliPAP 3-year warranty on humidification heater
Automatic Altitude Adjustment	Yes	DeVilbiss SmartLink® Compatible	Yes
DC Operable	Yes	SmartCode Feature	Yes
ABOB Analog Output for PSG	Yes (pressure only)		

* Tested per ISO 17510

** U.S. Patent 5,503,146



DeVilbiss Healthcare LLC | 100 DeVilbiss Drive | Somerset, PA 15501
800.338.1988 | 814.443.4881 | www.DeVilbissHealthcare.com