

PulmoNeb® Compressor Nebulizer System

Industry-leading
performance and reliability



DeVilbiss PulmoNeb
Compressor Nebulizer System
Model 3655X Series



The **DeVilbiss PulmoNeb Compressor Nebulizer** joins the top-selling compressor line in the world, building upon the **field-proven advantages** of the DeVilbiss Pulmo-Aide® Compact Compressor.

The PulmoNeb System offers **best-in-class performance** in a design that makes it **easy for patients** to operate. Plus, the unit is free of internal tubing to **improve reliability** and reduce the chance of leaks.

In addition, integral wire guides molded into the system's housing prevent wires from coming into contact with moving parts. Ball bearing connecting rod design further **enhances system reliability**.

Performance Features

- Maximum compressor pressure – **35 psig or greater**
- Compressor free air flow – **9 lpm or greater**
- Tested to **at least 2000 hours** of life
- **5-year** limited warranty
- **Easy-to-access**, top-of-unit tubing connector facilitates attachment/removal – especially important for patients with limited dexterity
- Easily accessible one-piece MicroPore® Filter **requires no tools** to change

Design Features

- **Easy to clean** and maintain
- Darker body color **maintains its appearance** after a prolonged period of use
- **Exclusive** simple cord-wrap design
- Metal frame motor construction
- Ball bearing connecting rod design
- Integral internal wire guides
- Robust tubing connector
- Peel-off serial number labels on the outside of each single carton
- **User video included**



Model 3655DX includes: VixOne Disposable Nebulizer, tubing, instruction guide, and patient instructional DVD.



Model 3655TX includes: VixOne Reusable Nebulizer, tubing, instruction guide, and patient instructional DVD. Also includes a disposable nebulizer.

VixOne™ Nebulizer

- MMAD of 2.7 microns
- Respirable Fraction – 84% below 5 microns
- At least 18% more medication delivered
- 7' kink-resistant flexhose tubing
- Latex-free
- Reusable model is dishwasher safe to promote frequent cleaning
- Low residual volume (0.7cc)
- High-capacity medication bowl (>10cc)
- Comfortable mouthpiece design
- Anti-spill free design provides secure seal
- Nebulizes in upright position or at angles up to 45° to enhance comfort during treatment
- Universal connector accommodates standard adult and pediatric masks, as well as the PediNeb™ pediatric pacifier mouthpiece

Accessories

3655D-601	Filter (5 pack)
3655D-602	Filter (25 pack)
3655PN-612	PediNeb Pacifier (12 pack)
900D-610	VixOne Reusable Nebulizer (case of 10)
3655D-621	VixOne Disposable Nebulizer (case of 50)
AAM12	Adult Aerosol Mask (12 pack)
AAM20	Adult Aerosol Mask (50 pack)
PAM10	Pediatric Aerosol Mask (50 pack)
DL1050	Pediatric Character Mask (50 pack)

Technical Specifications

Product Dimensions	7.5" H x 4.0" W x 6.8" D	Filter Replacement	No tools
Product Weight	4.1 lbs.	Warranty	5 years
Max. Air Pressure	35+ psig	Safety Approvals	ETL
Max. Liter Flow	9+ lpm	Motor/Pump Assy.	Primarily metal parts
MMAD	5 microns or less	Neb & Cord Storage	Yes
Sound	56 dBA	Master Carton Qty.	6



DeVilbiss Healthcare LLC | 100 DeVilbiss Drive | Somerset, PA 15501 USA
800.338.1988 | 814.443.4881 | www.DeVilbissHealthcare.com