

## Nasal-Aire<sup>®</sup> II CPAP Interface

### Something New is in the Aire...

The Nasal-Aire II is the **next generation of the ever-popular Nasal-Aire** interface and shares many of its **comfort and compliance benefits**, adding improved performance and a reduction in sound level.

This **lightweight interface is comfortable for patients, simple to use and easy to care for**. Nasal-Aire II is **sold with all six sizes for easy fitting** and includes headgear for stability.



Nasal-Aire II  
CPAP Interface  
(K2A)

### Performance Features

- Ergonomically designed to fit any face with virtually no leaks
- Virtually eliminates uncomfortable pressure points on the face allowing the interface to be worn for any duration without discomfort
- Secure fit design **allows the user to sleep in any position**
- Minimal design **reduces feelings of claustrophobia**
- New ultra-soft nasal prong material for improved patient comfort

## Design Features

- 100% latex free
- Medium or large size fits 90% of patients
- Can be worn over the ears or down the chest
- Small and quiet with virtually no noise
- No nasal bridge or forehead contact **eliminates nasal bridge sores and air leaks near eyes**
- Swivel tubing coupler to accommodate free movement and optimal performance during the night
- Supplemental oxygen ports

## Reimbursement Information

Product Description	Code Description	HCPCS Coding	Frequency	Allowable Range
Cannula	Cushion for nasal application device	A7032	<b>Billing at initial setup</b> Two per month	\$34 - \$40
Nasal-Aire II All Size, New Patient and Frame Kits	Nasal interface (mask or cannula type) used with positive airway pressure device, with or without head strap	A7034	<b>Billing at initial setup</b> One per 3 months	\$99 - \$117
Headgear	Headgear used with positive airway pressure device	A7035	<b>Billing at initial setup</b> One per 6 months	\$33 - \$39

## Accessories and Replacement Parts

PART NUMBER	PRODUCT DESCRIPTION	PART NUMBER	PRODUCT DESCRIPTION
K2A	Nasal-Aire II All Size Kit, Includes XL, LG, MP, MD, SM, XS*	SM02	Nasal-Aire II Cannula Only, SM
K2LG	Nasal-Aire II New Patient Setup Kit, MP, LG, XL*	XL06	Nasal-Aire II Cannula Only, XL
K2MD	Nasal-Aire II New Patient Setup Kit, SM, MD, MP*	XS01	Nasal-Aire II Cannula Only, XS
K2MP	Nasal-Aire II New Patient Setup Kit, MD, MP, LG*	X2105	Nasal-Aire II Headgear with OHS Assembly
K2SM	Nasal-Aire II New Patient Setup Kit, XSM, SM, MD*	P430	Nasal-Aire II Petite – Headgear
PA401	Nasal-Aire II Petite, Size A**	PA411	Nasal-Aire II Petite – Cannula Only, Size A
PB402	Nasal-Aire II Petite, Size B**	PB412	Nasal-Aire II Petite – Cannula Only, Size B
PC403	Nasal-Aire II Petite, Size C**	PC413	Nasal-Aire II Petite – Cannula Only, Size C
PD404	Nasal-Aire II Petite, Size D**	PD414	Nasal-Aire II Petite – Cannula Only, Size D
PE405A	Nasal-Aire II Petite, Size E**	PE415A	Nasal-Aire II Petite – Cannula Only, Size E
LG305	Nasal-Aire II, LG*	NAT320	Nasal-Aire II Frame Kit
MD303	Nasal-Aire II, MD*	NAT321	Nasal-Aire II Frame System, XS***
MP304	Nasal-Aire II, MP*	NAT322	Nasal-Aire II Frame System, SM***
SM302	Nasal-Aire II, SM*	NAT323	Nasal-Aire II Frame System, MD***
XL306A	Nasal-Aire II, XL*	NAT324	Nasal-Aire II Frame System, MP***
XS301	Nasal-Aire II, XS*	NAT325	Nasal-Aire II Frame System, LG***
LG05	Nasal-Aire II Cannula Only, LG	NAT326	Nasal-Aire II Frame System, XL***
MD03	Nasal-Aire II Cannula Only, MD	X2007	Nasal-Aire II Sizing Gauge
MP04	Nasal-Aire II Cannula Only, MP	X4007	Nasal-Aire II Petite Sizing Gauge

\* Includes headgear, frame, bolo and multiple nasal cannulas

\*\* Includes headgear, frame, nasal cannula and bolo

\*\*\* Includes cannula and tubing – no headgear



DeVilbiss Healthcare LLC | 100 DeVilbiss Drive | Somerset, PA 15501 USA  
800.338.1988 | 814.443.4881 | www.DeVilbissHealthcare.com