

HomeCare Suction Unit

A Fusion of Quality and Innovation



HomeCare Suction Unit
(Model 7305P-D)

For more than twenty years, DeVilbiss has been delivering the **best portable suction options**. The DeVilbiss 7305 Series Homecare Suction Unit **meets suction and capacity requirements in a small, lightweight design**.

The DeVilbiss 7305 series Homecare Suction Unit is a compact medical suctioning device which has been **designed for reliable, portable operation**. Because of the **small size, light weight and DC operation**, the DeVilbiss 7305 Suction Unit is ideal for providing suction in the home, or on the go with the **DC cord**. The 7305P series with **internal rechargeable battery** is also available. **Two container options** give the choice between the **800 ml disposable container** and optional long-term **reusable 1,200 ml container**.

Performance Features

- High Performance – Vacuum adjustments allow for 80-550 mm Hg and a free flow of 27 LPM
- Adjustable Flow Regulator – Easy-to-turn knob conveniently placed directly below gauge for easy use and control. Knob has safety lock feature which will not allow the knob to be twisted off.
- The Homecare Suction Unit 7305P series meets the RTCA/DO160D Aircraft Standard as well as the ISO-10079.1:1999 Suction Standard
- Increased Battery Capabilities – For use up to one hour. Provides the patient with added security when needed most
- The high flow capabilities (27 LPM) of the HomeCare Suction Unit can return a patient to a state of comfort quickly.

Design Features

- Easy-to-Read Gauge – The unit has conveniently placed gauge for easy reading. This is especially important when precise measurements of vacuum are needed for the pediatric patient.
- The 7305 series is equipped with a switch mode power supply allowing operation on any AC voltage (100-240 VAC, 50/60 Hz).
- Two-year warranty on Suction Unit
- Compact Design – All of the high-performance capabilities are conveniently placed in a small compact package weighing less than 6.3 lbs. The compact design with the high-performance

capabilities makes the 7305 series suction unit the smallest and strongest portable suction unit of its kind.

- Low battery and charge indicator lights on portable unit (7305P-D)

External Filter Configuration



- Connection tubing
- Filter
- Elbow
- Lid with float shutoff
- Jar
- 6' patient tubing

Internal Filter Configuration



- Connection tubing
- Elbow with integrated filter
- Float ball
- Lid
- Splash guard
- Jar
- 6' patient tubing

Available Configurations

PART NUMBER	PRODUCT DESCRIPTION
7305D-D	Suction Unit, 800cc disposable container with internal filter, tubing, elbow, and AC & DC power cords
7305D-D-EXF	Suction Unit, 800cc disposable container with external filter, tubing, elbow, and AC & DC power cords
7305P-D	Suction Unit, 800cc disposable container with internal filter, tubing, elbow, AC & DC power cords, and rechargeable battery
7305P-D-EXF	Suction Unit, 800cc disposable container with external filter, tubing, elbow, AC & DC power cords, and rechargeable battery

Accessories and Replacement Parts

PART NUMBER	PRODUCT DESCRIPTION
7305D-606	Carrying Case
7304D-619	12V DC power cord [1 each]
7305P-613	AC to DC adapter/charger
099HD-614	Hospital-grade power cord (120 VAC)
6710D-609	Power cord for US
7305P-631	Power cord for Continental Europe

Container Kits and Parts

Disposable Container Kits and Parts

PART NUMBER	PRODUCT DESCRIPTION
7305D-603	External Filter Disposable Container Complete Kit (800 ml disposable container, bacteria filter, elbow, 3-3/8" and 6' tubing)
7305D-602	External Filter Disposable Container (48 pack) (container and lid)
7305D-633	Collection Container Kit (800 ml disposable container, filter cartridge, 4-3/8" and 6' tubing package)
7305D-632	800 ml disposable container with filter cartridge and 4-3/8" tubing (48 each)
7305D-635	Filter cartridge (12 pack) (For Disposable Container)
7305D-608	External filter (12 pack)
6305D-611	6' Patient Tubing (container and lid)

Reusable Container Kits and Parts

PART NUMBER	PRODUCT DESCRIPTION
7314D-603	Collection Container Kit (1200 ml reusable container, bacteria filter, elbow, 3-3/8" tubing)
7314D-604	1200 ml reusable container (bacteria filter, elbow, 3-3/8" tubing) [6 pack]
7305D-608	External filter (12 pack)

Technical Specifications

Dimensions	9" x 7.0" x 8.0" [22.9cm x 17.8cm x 20.3cm]
7305D Series	3.8 lbs (1.7 kg)
7305P Series	6.3 lbs (2.9 kg)
Electrical Requirements	100-240 V AC 50/60 Hz .75 A max; 12 V DC, 33 W max
7305P Series Internal Rechargeable Battery	Factory Equipped
7305D Series Internal Rechargeable Battery	Not Factory Equipped (provider-installed option)
7305P & 7305D Series Vacuum Range	80 to 550 mm Hg
Air Flow @ pump inlet:	27 LPM (free flow) typical (may be less when running from internal battery)
7305 D&P Series Container Capacity	800 ml (cc) Disposable
7305 D&P Series Container Bottle Capacity-Optional	1,200 ml (cc) Reusable
7305P Series Warranty	Two-years limited, excluding internal battery and container
7305D Series Warranty	Two-years limited, excluding container
Internal Battery Warranty	90-day
7305D & 7305P Series Approvals	IEC 601-1; CAN/CSA-C22.2 No. 601.1-M90; UL 60601-1, EN 60601-1-2
7305P meets RTCA/DO & 160D Approvals for battery operation only; Airline use.	DO-160D - section 21 Category M, ISO-10079.1:1999



DeVilbiss Healthcare LLC | 100 DeVilbiss Drive | Somerset, PA 15501 USA
800.338.1988 | 814.443.4881 | www.DeVilbissHealthcare.com