

## Nebulizers, Masks and Accessories

### Efficient solutions for delivering aerosol therapy

The DeVilbiss compressor accessory line offers several products designed to provide efficient solutions to common issues. For the pediatric patient, the Pediatric Character Mask ensures that children receive proper treatment.

This child-friendly mask features a swivel connector for aerosol treatment in almost every position, while its custom vent system keeps aerosol out of the child's eyes. For infants, the patented PediNeb™ Pacifier system, developed by a nurse practitioner, provides nasal delivery of aerosol while the pacifier soothes the patient.

The VixOne nebulizer, available in reusable and disposable models, is the economical, high-performance choice in small-volume nebulizers. Its ideal particle size production allows for high-volume lung deposition, which results in maximum therapeutic effect by increasing the amount of medication delivered deep within the lungs. The corrugated holding chamber, standard with every unit, means patients receive 18 percent more aerosol than with a traditional open system.



Pediatric  
Character Mask



PediNeb  
Pacifier



VixOne Reusable  
Nebulizer

**CE** 0044  
These products are latex-free

### DeVilbiss Breathies Bunch™

This unique educational program is packed with important information, including a guide to childhood asthma, log book, entertaining videos, colorful stickers and more. The Breathies Bunch Website also complements the materials and includes interactive games to educate children about asthma and keep them compliant with their therapy.

## VixOne™ Nebulizer

- MMAD of 2.7 microns
- Respirable Fraction – 84% below 5 microns
- At least 18% more medication delivered
- 7' kink-resistant flexhose tubing
- Latex-free
- Reusable model is dishwasher safe to promote frequent cleaning
- Low residual volume (0.7cc)
- High-capacity medication bowl (>10cc)
- Comfortable mouthpiece design
- Anti-spill free design provides secure seal
- Nebulizes in upright position or at angles up to 45° to enhance comfort during treatment
- Universal connector accommodates standard adult and pediatric masks, as well as the PediNeb™ pediatric pacifier mouthpiece



## PediNeb Pacifier

The PediNeb was specially designed to enhance the delivery of aerosolized medication to infants and toddlers while reducing the apprehension experienced by these young patients.

- Pacifier calms the child during aerosol treatment
- Medication port directs aerosolized medication to the nose and reduces medication waste associated with the “blow-by” technique
- Easy to operate, clean and maintain
- Patented design was developed by a pediatric nurse practitioner

## Ordering Accessories

ITEM	DEVILBISS PART #	QUANTITY
VixOne Disposable Nebulizer	3655D-621	50 per case
VixOne Reusable Nebulizer	900D-610	10 per case
PediNeb Pacifier	3655PN-612	12 per case
Pediatric Character Mask	DL1050	50 per case
Adult Aerosol Mask	AAM12	12 per case
	AAM20	50 per case
Pediatric Aerosol Mask	PAM10	50 per case
Peak Flow Meter	600-112	12 per case
Breathies Bunch Compliance Kit	DBB-kids	Each

## Aerosol Masks



- Soft, pliable, latex-free masks conform comfortably around a wide array of facial contours.



DeVilbiss Healthcare LLC | 100 DeVilbiss Drive | Somerset, PA 15501 USA  
 800.338.1988 | 814.443.4881 | www.DeVilbissHealthcare.com