

DeVilbiss IntelliPAP® Specifications

Specification	Standard	Standard Plus	AutoAdjust®	Bilevel	AutoBilevel
Product Number	DV51D	DV53D	DV54D	DV55D	DV57D
Product Dimensions	6.4" H x 6.5" W x 8.4" D	6.4" H x 6.5" W x 8.4" D	6.4" H x 6.5" W x 8.4" D	6.4" H x 6.5" W x 8.4" D	6.4" H x 6.5" W x 8.4" D
Product Weight	4.45 lb with humidifier and PAP	4.45 lb with humidifier and PAP	4.45 lb with humidifier and PAP	4.45 lb with humidifier and PAP	4.45 lb with humidifier and PAP
Pressure Range	3-20 cm H ₂ O (.5 increments)	3-20 cm H ₂ O (.5 increments)	3-20 cm H ₂ O (.5 increments)	3-25 cm H ₂ O (.5 increments)	3-25 cm H ₂ O (.5 increments)
Universal Power Supply	Yes	Yes	Yes	Yes	Yes
AC Cord Length	8 feet	8 feet	8 feet	8 feet	8 feet
Automatic Altitude Adjustment	Yes	Yes	Yes	Yes	Yes
DC Operable	Yes	Yes	Yes	Yes	Yes
ABOB Analog Output For PSG	Yes (pressure only)	Yes	Yes	Yes	Yes
Automatic Leak Compensation	Yes	Yes	Yes	Yes	Yes
Auto On/Off	Yes (can be enabled or disabled)	Yes (can be enabled or disabled)	Yes (can be enabled or disabled)	Yes (can be enabled or disabled)	Yes (can be enabled or disabled)
Motor Hour Meter	Yes	Yes	Yes	Yes	Yes
Therapy Compliance Meter	Yes	Yes	Yes	Yes	Yes
Standard Warranty	3-year warranty on IntelliPAP	3-year warranty on IntelliPAP	3-year warranty on IntelliPAP	3-year warranty on IntelliPAP	3-year warranty on IntelliPAP
DeVilbiss SmartLink® Compatible	Yes	Yes	Yes	Yes	Yes
SmartCode® Feature	Yes	Yes	Yes	Yes	Yes
Tubing Lengths	6 feet or 10 feet	6 feet or 10 feet	6 feet or 10 feet	6 feet or 10 feet	6 feet or 10 feet



DeVilbiss Healthcare | 100 DeVilbiss Drive | Somerset, PA 15501 | USA
800.338.1988 | 814.443.4881 | www.DeVilbissHealthcare.com