

Convenient.

For more than half a century, DeVilbiss Healthcare has provided unmatched safety and performance in the airway suction device category, and we're confident this is one of our best units yet. In addition to being quiet, portable and powerful, the sleek DeVilbiss Vacu-Aide QSU includes these convenient features:

- ◀ **Two-year warranty on suction unit**
- ◀ **Low battery and charge indicator lights to help prevent interruptions in use**
- ◀ **Conveniently placed flow regulator knob for easy use and control**
- ◀ **Innovative container assembly with integrated bacteria filter**
- ◀ **Container and filter design to maintain appropriate flow rate to meet ISO guidelines**
- ◀ **Float shutoff on container to protect against overflow**
- ◀ **Integrated bacteria filter to prevent back flow**

Available Configurations

PART NUMBER	PRODUCT DESCRIPTION
7314P-D	Vacu-Aide QSU, 800 ml disposable container, internal battery, carrying case, 6' tubing, AC/DC adapter and power cord.
7314D-D	Vacu-Aide QSU, 800 ml disposable container, 6' tubing, AC/DC adapter and power cord. No internal battery.
7314P-D-EXF	Vacu-Aide QSU, 800 ml disposable container with external filter, internal battery, 6' tubing, AC/DC adapter and 12V power cord.
7314D-D-EXF	Vacu-Aide QSU, 800 ml disposable container with external filter, 6' tubing, AC/DC adapter and 12V power cord. No internal battery.



DeVilbiss Healthcare LLC | 100 DeVilbiss Drive | Somerset, PA 15501 USA
800.338.1988 | 814.443.4881 | www.DevilbissHealthcare.com

© 2014 DeVilbiss Healthcare LLC 07.14. All Rights Reserved.
DeVilbiss® and Vacu-Aide® are registered trademarks of DeVilbiss Healthcare. LT-2180 Rev D

Vacu-Aide® QSU Quiet Suction Unit



Quiet. Portable. Powerful.



Quiet.

Without the added noise that can often intensify feelings of discomfort, the DeVilbiss Vacu-Aide QSU quickly, efficiently and safely evacuates secretions from a patient's airway. In fact, with a 50% reduction in sound, it is the quietest high flow/high suction portable unit on the market.

Many times, suctioning in public is necessary, especially while traveling. For that reason, having a suction unit with a lower sound level is a must. In addition, when your situation demands a quiet environment, there is no question that the discreet Vacu-Aide QSU can help keep it that way.



Vacu-Aide QSU (7314P-D) includes the unit, internal battery, carrying case, 6' tubing, AC/DC adapters and power cord. Also available without internal battery (7314D-D).

Portable.

When a patient needs a suction unit to clear airway secretions, suctioning can be required anytime and anywhere. The portable and lightweight Vacu-Aide QSU is easily transported and can accommodate travel to any country. With a multitude of power options, features and accessories, it can be used virtually anywhere under some of the most extreme circumstances.

- Small footprint and lightweight design – only 6.6 lbs.
- Discreet carrying case can be used during operation
- Optional sealed lead acid battery for use up to one-hour operation
- Switch mode power supply, allowing operation on any AC voltage (100-240 VAC, 50/60 Hz)

50% Quieter

Powerful.

DeVilbiss has long been known for designing and manufacturing quality, durable respiratory medical equipment, and the Vacu-Aide QSU joined the ranks as a high performance unit. This powerful unit with adjustable vacuum range makes it a smart choice for patients from newborn to adult.

- Vacuum adjustment allows for 50-550 mm Hg
- Free flow of 27 liters/minute
- Meets ISO 10079-1:2009 high vacuum, high flow classification and can return a patient to a state of comfort quickly
- Meets American Association for Respiratory Care (AARC) guidelines for neonatal, infant, child and adult homecare suctioning

Quiet. Portable. Powerful.

To learn more, visit
www.DeVilbissHealthcare.com