

## Nasal-Aire<sup>®</sup> II CPAP Interface

### Something New is in the Aire...

The Nasal-Aire II is the **next generation of the ever-popular Nasal-Aire** interface. The Nasal-Aire II shares many of the **comfort and compliance benefits** of the original Nasal-Aire, adding improved performance and a reduction in sound level.

**This lightweight interface is comfortable for patients, simple to use and easy to care for.** Nasal-Aire II is **sold with all six sizes for easy fitting** and includes headgear for stability.



Nasal-Aire II  
CPAP Interface  
(K2A)

### Performance Features

- Ergonomically designed to fit any face with virtually no leaks
- Virtually eliminates uncomfortable pressure points on the face allowing the interface to be worn for any duration without discomfort
- Secure fit design **allows the user to sleep in any position**
- Minimal design **reduces feelings of claustrophobia**
- New ultra-soft nasal prong material for improved patient comfort

## Design Features

- 100% latex free
- Medium or large size fits ninety percent of patients
- Can be worn over the ears or down the chest
- Small and quiet with virtually no noise
- No nasal bridge or forehead contact **eliminates nasal bridge sores and air leaks near eyes**
- Swivel tubing coupler to accommodate free movement and optimal performance during the night
- Supplemental oxygen ports

## Reimbursement Information

| Product Description       | Code Description   | HCPCS Coding | Frequency   | Allowable Range |
|---------------------------|--|--------------|---|-----------------|
| NASAL-AIRE II INSERT      | Cushion for nasal application device   | A7032        | <b>Billing at initial setup</b><br>Two per month    | \$34 -\$40      |
| NASAL-AIRE II APPLICATION | Nasal interface (mask or cannula type) used with positive airway pressure device, with or without head strap | A7034        | <b>Billing at initial setup</b><br>One per 3 months | \$99 -\$117     |
| HEADGEAR                  | Headgear used with positive airway pressure device   | A7035        | <b>Billing at initial setup</b><br>One per 6 months | \$33 -\$39      |

## Accessories and Replacement Parts

| PART NUMBER | PRODUCT DESCRIPTION  | PART NUMBER | PRODUCT DESCRIPTION                         |
|-------------|--|-------------|---|
| K2A         | Nasal-Aire II All Size Kit, Includes XL, LG, MP, MD, SM, XS* | LG05        | Nasal-Aire II Cannula Only, LG              |
| K2LG        | Nasal-Aire II New Patient Setup Kit, MP, LG, XL*             | MD03        | Nasal-Aire II Cannula Only, MD              |
| K2MD        | Nasal-Aire II New Patient Setup Kit, SM, MD, MP*             | MP04        | Nasal-Aire II Cannula Only, MP              |
| K2MP        | Nasal-Aire II New Patient Setup Kit, MD, MP, LG*             | SM02        | Nasal-Aire II Cannula Only, SM              |
| K2SM        | Nasal-Aire II New Patient Setup Kit, XSM, SM, MD*            | XL06        | Nasal-Aire II Cannula Only, XL              |
| PA401       | Nasal-Aire II Petite, Size A                                 | XS01        | Nasal-Aire II Cannula Only, XS              |
| PB402       | Nasal-Aire II Petite, Size B                                 | NAT320      | Nasal-Aire II Frame Kit                     |
| PC403       | Nasal-Aire II Petite, Size C                                 | X2105       | Nasal-Aire II Headgear with OHS             |
| PD404       | Nasal-Aire II Petite, Size D                                 | Assembly    |   |
| PE405A      | Nasal-Aire II Petite, Size E                                 | P430        | Nasal-Aire II Petite – Headgear             |
| LG305       | Nasal-Aire II, LG*   | PA411       | Nasal-Aire II Petite – Cannula Only, Size A |
| MD303       | Nasal-Aire II, MD*   | PB412       | Nasal-Aire II Petite – Cannula Only, Size B |
| MP304       | Nasal-Aire II, MP*   | PC413       | Nasal-Aire II Petite – Cannula Only, Size C |
| SM302       | Nasal-Aire II, SM*   | PD414       | Nasal-Aire II Petite – Cannula Only, Size D |
| XL306A      | Nasal-Aire II, XL*   | PE415A      | Nasal-Aire II Petite – Cannula Only, Size E |
| XS301       | Nasal-Aire II, XS*   | X2007       | Sizing Gauge – Nasal-Aire II                |
|             |  | X4007       | Sizing Gauge – Nasal-Aire II Petite         |

\* Includes headgear, frame, bolo and multiple nasal cannulas.

\*\* Includes headgear, frame, nasal cannula and bolo.



DeVilbiss Healthcare LLC | 100 DeVilbiss Drive | Somerset, PA 15501 USA  
800.338.1988 | 814.443.4881 | [www.DeVilbissHealthcare.com](http://www.DeVilbissHealthcare.com)