



DeVilbiss PulmoNeb™ Compressor Nebulizer System

Industry-leading performance and reliability

The new **DeVilbiss PulmoNeb Compressor Nebulizer** joins the top-selling compressor line in the world, building upon the **field-proven advantages** of the DeVilbiss Pulmo-Aide® Compact Compressor.

The PulmoNeb System (including a children's version) offers **best-in-class performance** in a highly intuitive design that makes it **easy for patients** to operate. Plus, the unit is free of internal tubing to **improve reliability** and reduce the chance of leaks.

In addition, integral wire guides molded into the system's housing prevent wires from coming into contact with moving parts. Ball bearing connecting rod design further **enhances system reliability**.



DeVilbiss PulmoNeb
Compressor Nebulizer Unit
Model 3655TX



Performance Features

- Maximum compressor pressure – **35 psig or greater**
- Compressor free air flow – **9 lpm or greater**
- Tested to **at least 2000 hours** of life
- **5-year** limited warranty
- **Easy-to-access**, top-of-unit tubing connector facilitates attachment/removal – especially important for patients with limited dexterity
- Easily accessible one-piece MicroPore® Filter **requires no tools** to change

Design Features

- **Easy to clean** and maintain
- Darker body color **maintains its appearance** after a prolonged period of use
- **Exclusive** simple cord-wrap design
- Metal frame motor construction
- Ball bearing connecting rod design
- Integral internal wire guides
- Robust tubing connector
- Peel-off serial number labels on the outside of each single carton
- **User video included**

800 Reusable Jet Nebulizer

- MMAD 3.6-4 microns
- Approximately 75% of particles are below 5 microns in size
- Reduces treatment time nearly in half
 - 3cc treatment in 7-8 minutes
- 6 cc capacity
- Dishwasher safe to promote frequent cleaning
- Ergonomic design fits comfortably in the hand
 - even with limited functionality
- 4 components for ease of assembly/disassembly
- Comfortable mouthpiece design
- Standard 22mm fitting, no adapters required
- Can be used inline with a ventilator
- 1 year warranty
- Tested for use with Albuterol, Pulmicort (budesonide), Ipratropium bromide.
- Suitable for use with most common liquid medications



Model 3655DX includes: Disposable nebulizer, tubing, instruction guide, and patient instructional DVD.



Model 3655TX includes: DeVilbiss 800 Reusable Jet nebulizer, tubing, instruction guide, and patient instructional DVD. Also includes a disposable nebulizer.



Model 3655PL-KIDS includes: Disposable nebulizer, tubing, instruction guide, and Breathies Bunch program materials including resource guide, log book, stickers, Breathies nebulizer hugger and instructional DVD. Refer to sell sheet LT-XXXX for more info. Components may vary.

Accessories

3655D-601	Filter (5 pack)
3655D-602	Filter (25 pack)
800D	Reusable Nebulizer – blue (1 each)
800C	Reusable Nebulizer – purple/green (1 each)
4650D-621	Disposable Nebulizer (50 pack)
AAM12	Adult Aerosol Mask (12 pack)
AAM20	Adult Aerosol Mask (50 pack)
PAM10	Pediatric Aerosol Mask (50 pack)
DLI050	Pediatric Children's Character Mask (50 pack)

Technical Specifications

Product Dimensions	7.5" H x 4.0" W x 6.8" D	Filter Replacement	No tools
Product Weight	4.1 lbs.	Warranty	5 years
Max. Air Pressure	35+ psig	Safety Approvals	ETL
Max. Liter Flow	9+ lpm	Motor/Pump Assy.	Primarily metal parts
MMAD	5 microns or less	Neb & Cord Storage	Yes
Sound	56 dBA	Master Carton Qty.	6



DeVilbiss Healthcare | 100 DeVilbiss Drive | Somerset, PA 15501
800.338.1988 | 814.443.4881 | www.DeVilbissHealthcare.com