

Frequently Asked Questions VacuAide[®] Compact Suction Unit

Q1: How long will the battery operate the unit?

A1: Up to one-hour

Q2: What is recommended for the initial charge of the battery?

A2: Fully charge the battery for a minimum of five hours

Q3: If not in use for an extended period of time are there any recommendations for keeping the battery in good condition?

A3: Yes, it's our recommendation that the battery be recharged a minimum of every 6 months

Q4: If traveling internationally is there anything needed to use my Compact Suction Unit?

A4: Yes, the correct power cord must be used, as the unit can operate on any AC voltage.

Q5: How do I get replacement accessories and parts for my Compact Suction Unit?

A5: Through an Authorized DeVilbiss Provider, which can be found on our web site under Find a Provider, located at the very top of our web pages.

Q6: What is the warranty of my Compact Suction Unit?

A6: 7310 series – 2 years from date of purchase (excludes the battery & collection canister); the internal rechargeable battery is warranted for 6 months.

