

## Frequently Asked Questions VacuAide<sup>®</sup> Compact Suction Unit

**Q1: How long will the battery operate the unit?**

*A1: Up to one-hour*

**Q2: What is recommended for the initial charge of the battery?**

*A2: Fully charge the battery for a minimum of five hours*

**Q3: If not in use for an extended period of time are there any recommendations for keeping the battery in good condition?**

*A3: Yes, it's our recommendation that the battery be recharged a minimum of every 6 months*

**Q4: If traveling internationally is there anything needed to use my Compact Suction Unit?**

*A4: Yes, the correct power cord must be used, as the unit can operate on any AC voltage.*

**Q5: How do I get replacement accessories and parts for my Compact Suction Unit?**

*A5: Through an Authorized DeVilbiss Provider, which can be found on our web site under Find a Provider, located at the very top of our web pages.*

**Q6: What is the warranty of my Compact Suction Unit?**

*A6: 7310 series – 2 years from date of purchase (excludes the battery & collection canister); the internal rechargeable battery is warranted for 6 months.*

