

## Built to last. Built for life.

When you're on the go, and need supplemental oxygen, you can go as you please with the DeVilbiss iGo Portable Oxygen System – the rugged and reliable choice in portable oxygen systems.

The iGo System has been performance-tested to meet the demands of your active lifestyle. In fact, no other system of its kind is built to withstand the punishment that iGo takes in stride. The compact and lightweight iGo System gives you the freedom to expand your travel horizons.



## DeVilbiss Is With You Every Step of the Way

To help you acclimate to your new iGo Portable Oxygen System, DeVilbiss offers a highly informative DVD that helps you get the most out of your oxygen therapy, while demonstrating how to properly use and care for your system.

To learn more,  
call 800.338.1988 or visit  
[www.iGopoc.com](http://www.iGopoc.com)



DeVilbiss Healthcare LLC | 100 DeVilbiss Drive | Somerset, PA 15501 USA  
800.338.1988 | 814.443.4881 | [www.DeVilbissHealthcare.com](http://www.DeVilbissHealthcare.com)

## DeVilbiss iGo® Portable Oxygen System



FAA  
Approved  
for air travel!



## Choice of oxygen delivery modes.

Exclusive PulseDose® oxygen-conserving technology from DeVilbiss delivers brief and consistent bursts of oxygen – even at higher breath rates. For many, these short bursts are almost undetectable and more comfortable than continuous operation. Plus, PulseDose helps reduce throat and nasal dryness, and extends battery life too.

The iGo System also offers a continuous flow mode that is recommended for use during sleep.

## Dependable, fast-charging lithium batteries.

For added travel flexibility, the iGo battery quickly recharges in 2 to 4.5 hours and automatically recharges when connected to an AC power source. If the power source is ever interrupted, your iGo System automatically switches to battery operation.

## A new standard of accuracy and performance.

Every iGo System is factory-tested for proper oxygen purity, flow rate, pulse volume and battery run time to let you travel with confidence – now and in the future. Extensive independent testing has proven the unit's effectiveness.

## Safety first.

The iGo System's patented oxygen sensing device helps ensure accurate oxygen delivery and reduces periodic maintenance. In addition, audible beeps alert you to oxygen output and power problems, including low battery plus any unit malfunction.

Rugged rolling carry case with wheels and pull-up handle easily navigates bumps and rough pavement

Sturdy construction provides years of problem-free operation

Connects to 50-foot tubing/cannula in continuous flow mode; 35-foot tubing in PulseDose mode

Weighs less than 20 pounds with battery installed

Settings 1 to 6 in PulseDose mode; 1 to 3 lpm in Continuous Flow mode

Up to 5.4 hours of extended battery life in PulseDose mode

System includes option of 1 or 2 batteries, plus home AC power cord and car DC power cord

Full 3-year limited warranty (1-year limited warranty on battery)

Easy-access filter compartment

Easy-to-read battery status gauge

Control panel can be easily viewed – even under an airplane seat



The DeVilbiss iGo fits in your car for easy traveling.

